BORIS promotes chromatin regulatory interactions in treatment-resistant cancer cells

David N. Debruyne^{1,2,15}, Ruben Dries^{1,2,3,15}, Satyaki Sengupta^{1,2}, Davide Seruggia^{1,4,5,6}, Yang Gao^{1,2}, Bandana Sharma^{1,2}, Hao Huang^{1,2}, Lisa Moreau⁷, Michael McLane^{1,2}, Daniel S. Day^{8,9}, Eugenio Marco^{3,10}, Ting Chen¹¹, Nathanael S. Gray^{12,13}, Kwok-Kin Wong¹⁴, Stuart H. Orkin^{1,4,5,6}, Guo-Cheng Yuan^{3,10}, Richard A. Young^{8,9} & Rani E. George^{1,2,4}*

The CCCTC-binding factor (CTCF), which anchors DNA loops that organize the genome into structural domains, has a central role in gene control by facilitating or constraining interactions between genes and their regulatory elements^{1,2}. In cancer cells, the disruption of CTCF binding at specific loci by somatic mutation^{3,4} or DNA hypermethylation⁵ results in the loss of loop anchors and consequent activation of oncogenes. By contrast, the germ-cell-specific paralogue of *CTCF*, *BORIS* (brother of the regulator of imprinted sites, also known as *CTCFL*)⁶, is overexpressed in several cancers^{7–9}, but its contributions to the malignant phenotype remain unclear. Here we show that aberrant upregulation of *BORIS* promotes chromatin interactions in *ALK*-mutated, *MYCN*-amplified neuroblastoma¹⁰ cells that develop resistance to ALK inhibition. These cells are reprogrammed to a distinct phenotypic state during the acquisition of resistance, a process defined by the initial loss of *MYCN* expression followed by subsequent overexpression of *BORIS* and a concomitant switch in cellular dependence from MYCN to BORIS. The resultant BORIS-regulated alterations in chromatin looping lead to the formation of super-enhancers that drive the ectopic expression of a subset of proneural transcription factors that ultimately define the resistance phenotype. These results identify a previously unrecognized role of BORIS—to promote regulatory chromatin interactions that support specific cancer phenotypes.




Fig. 1 | Targeted therapy resistance in neuroblastoma is associated with transcriptional reprogramming and a switch in dependency from amplified MYCN to BORIS. a, Top, schematic representation of the development of resistance. Bottom, dose-response curves of TAE684-sensitive and -resistant Kelly neuroblastoma cells incubated in increasing concentrations of TAE684 for 72 h. Data are mean \pm s.d., n = 3 biological replicates. **b**, Heat map of gene expression values in sensitive versus resistant cells (n = 2 biological replicates). Rows are z-scores calculated for each gene in both cell types. c, PCA of scRNA-seq data of sensitive (n = 5,432), intermediate resistant (IR; n = 6,376) and resistant (n = 6.379) cells showing the first two principal components (PCs). d, Pseudotime analysis of transcription factor expression during the development of resistance. e, ChIP-seq signals of genome-wide MYCN binding in sensitive and resistant cells, reported as reads per million (RPM) per base pair (bp) for each chromosome (chr). f, PCA of gene expression profiles showing the first two principal components (n = 2biological replicates). g, Dose-response curves for TAE684 (half-maximum inhibitory concentration (IC₅₀) values in parenthesis) and immunoblot analysis (representative of two independent experiments) of BORIS and MYCN expression in sensitive cells expressing short hairpin RNA (shRNA) against MYCN (MYCN^{KD}) and doxycycline-inducible BORIS (BORIS^{Ind}), treated with dimethylsulfoxide (DMSO) or 1 µM TAE684, with or without doxycycline (DOX). Data are mean \pm s.d., n = 3 biological replicates.

¹Department of Pediatric Oncology, Dana-Farber Cancer Institute, Boston, MA, USA. ²Department of Pediatrics, Harvard Medical School, Boston, MA, USA. ³Department of Biostatistics and Computational Biology, Dana-Farber Cancer Institute, Boston, MA, USA. ⁴Division of Hematology/Oncology, Boston Children's Hospital, Boston, MA, USA. ⁵Harvard Stem Cell Institute, Harvard Medical School, Boston, MA, USA. ⁶Howard Hughes Medical Institute, Boston, MA, USA. ⁷Department of Radiation Oncology, Dana-Farber Cancer Institute, Boston, MA, USA. ⁸Whitehead Institute for Biomedical Research, Cambridge, MA, USA. ⁹MIT Department of Biology, Cambridge, MA, USA. ¹⁰Department of Biostatistics, Harvard TC Chan School of Public Health, Boston, MA, USA. ¹¹Department of Medical Oncology, Dana-Farber Cancer Institute, Boston, MA, USA. ¹²Department of Cancer Biology, Dana-Farber Cancer Institute, Boston, MA, USA. ¹¹Department of Medical Oncology, Dana-Farber Cancer Institute, Boston, MA, USA. ¹²Department of Biology, Cambridge, MA, USA. ¹³Department of Biological Chemistry and Molecular Pharmacology, Harvard Medical School, Boston, MA, USA. ¹⁴Division of Hematology and Medical Oncology, Laura and Isaac Perlmutter Cancer Center, New York University Langone Medical Center, New York, NY, USA. ¹⁵These authors contributed equally: David N. Debruyne, Ruben Dries. ^{*e-mail}: rani_george@dfci.harvard.edu

Unlike CTCF, which is uniformly expressed in healthy tissues and cancer cells, the expression of BORIS is typically restricted to the testis⁶ and embryonic stem cells¹¹ (Extended Data Fig. 1a). However, when aberrantly expressed in cancer⁷⁻⁹, it is associated with high-risk features that include resistance to treatment (Extended Data Fig. 1b, c). We identified BORIS as one of the most differentially expressed genes in neuroblastoma cells driven by amplified MYCN¹² and ALK(F1174L)¹³ and rendered resistant to ALK inhibition. Kelly human neuroblastoma cells were exposed to increasing concentrations of the ALK inhibitor TAE684¹⁴ until stable resistance was achieved (Fig. 1a, Extended Data Fig. 2a-d). The acquisition of stable resistance coincided not only with the loss of ALK phosphorylation-which indicates that the cells no longer required activation of this receptor tyrosine kinase to maintain their oncogenic properties-but also with the absence of other common instigators of resistance (Extended Data Fig. 2a, e-h; Supplementary Note 1). However, comparison of the gene expression profiles of the TAE684-sensitive and resistant cells showed generalized downregulation of transcription in the resistant cells, but with marked upregulation of a subset of transcription factors not typically associated with neuroblastoma cells^{15,16} (Fig. 1b).

We therefore proposed that the resistant cells had probably undergone transcriptional reprogramming during the development of resistance. To determine the dynamics of resistance development, we performed single-cell RNA sequencing (scRNA-seq) analysis on sensitive, intermediate and fully resistant cell states (Extended Data Fig. 3a). Principal component analysis (PCA) indicated a stepwise transition as cells progressed from the sensitive to the fully resistant state (Fig. 1c). This transition was confirmed by distributed stochastic neighbour embedding (t-SNE)¹⁷, which clustered the cells into three nonoverlapping categories (Extended Data Fig. 3b, c). Pseudotime analysis based on the transcription factors that were differentially expressed throughout the development of resistance revealed that the initial major alteration was loss of MYCN expression, which persisted in stably resistant cells (Fig. 1d, Extended Data Fig. 3d, e). To understand this unexpected result, we analysed the status of MYCN in these cells, and found that although genomic amplification was retained, the MYCN locus was epigenetically repressed (Extended Data Fig. 3f, g). This state was accompanied by a genome-wide reduction of MYCN binding to DNA and a consequent revision of associated downstream transcription outcomes^{15,18,19} (Fig. 1e, Extended Data Fig. 3h). Coincident with this loss of transcriptional activity, the resistant cells were no longer dependent on MYCN for survival, unlike their sensitive controls, which underwent apoptosis after depletion of MYCN (Extended Data Fig. 3i). Subsequent resistance stages were defined by a gradual increase in the expression of the neural developmental markers SOX2 and SOX9²⁰, followed by upregulation of BORIS, ultimately leading to a fully resistant state in which BORIS expression was highest and detectable in essentially all cells (Fig. 1d, Extended Data Fig. 3j, k). Overexpression of BORIS, which coincided with promoter hypomethylation (Extended Data Fig. 4a, b), was also observed in additional neuroblastoma cell lines rendered resistant to TAE684 (SK-N-SH) or the CDK12 inhibitor E9²¹ (SK-N-BE(2)) (Extended Data Fig. 4c, d), which suggests that our findings are not restricted to a single cell line or kinase inhibitor. Indeed, overexpression of BORIS in tumours was significantly associated with high-risk disease and a poor outcome in patients with neuroblastoma treated with a variety of regimens (Extended Data Fig. 4e-g).

To clarify the role of BORIS in the resistance phenotype, we depleted its expression in resistant cells, and observed a partial reversal to the sensitive-cell state with re-emergence of MYCN and ALK expression (Fig. 1f, Extended Data Fig. 5a–c). However, this outcome was insufficient to maintain cell growth, as depletion of BORIS in resistant cells ultimately decreased cell viability (Extended Data Fig. 5d, e), which indicates a switch from MYCN to BORIS dependency with stable resistance. This transition was associated with changes in cellular growth kinetics—from a highly proliferative, *MYCN*-overexpressing sensitive state to an intermediate, slow-cycling phenotype that was partially reversed in fully resistant cells, coincident with overexpression of




Fig. 2 | *BORIS* overexpression is associated with its increased chromatin occupancy in resistant cells, whereas CTCF binding is unchanged. a, Scatter plot of BORIS binding in sensitive (Sens) and resistant (Res) cells for all detected BORIS-binding sites. BORIS peaks unique to resistant cells (n = 21,805; 91%), sensitive cells (n = 1,125; 4.7%) and shared between the two cell types (n = 1,086; 4.5%) are shown. b, Scatter plot of CTCF binding in sensitive and resistant cells for all detected CTCF-binding sites. CTCF peaks unique to resistant cells (n = 6,808; 8.3%), sensitive cells (n = 19,129; 23.2%) and shared between the two cell types (n = 56,438; 68.5%) are shown. c, Overlap between BORIS peaks that are unique to resistant cells and CTCF peaks shared between resistant and sensitive cells (top), and between resistant cell-specific BORIS peaks and sensitive cell-specific CTCF peaks (bottom). d, Meta-analysis of average ChIP-seq signals at resistant cell-specific BORIS-binding sites. All panels, n = 2 biological replicates.

BORIS (Extended Data Fig. 5f–h). Given the many sequential steps involved in the evolution of resistance, overexpression of *BORIS* alone was not adequate to induce this phenotype (data not shown). Instead, concomitant downregulation of *MYCN* expression and *BORIS* overexpression in the presence of ALK inhibition were required to generate resistance in sensitive cells (Fig. 1g). This combination of factors also led to increased expression of the transcription factors that were upregulated in the original TAE684-resistant cells, including *SOX2* and *SOX9* (Extended Data Figs. 3d, 5i). Thus, resistance to inhibition of ALK in neuroblastoma cells evolves through a multistep process that promotes a dependency switch from a dominant oncogenic stimulus—amplified *MYCN*—to a phenotypically distinct state characterized by overexpression of *BORIS*. In this context, the resistant cells ultimately become dependent on BORIS for survival, which supports a key role for this protein in maintenance of the resistance state.

We next asked whether the aberrant expression of BORIS, a DNAbinding protein⁶, affected its genome-wide occupancy in resistant cells. We observed a large (tenfold) gain in BORIS-bound peaks after chromatin immunoprecipitation followed by high-throughput sequencing (ChIP-seq) analysis in resistant cells: 22,891 versus 2,211 in sensitive cells (Fig. 2a, Extended Data Fig. 6a, b). By contrast, CTCF binding did not change substantially between sensitive and resistant cells (75,567 versus 63,246 peaks) (Fig. 2b). A considerable portion (n = 17,042;78%) of the BORIS peaks unique to resistant cells overlapped with CTCF peaks shared by both cell types (Fig. 2c), consistent with their heterodimerization²² (Extended Data Fig. 6c). However, only a small proportion (n = 1,903; 8.7%) overlapped with CTCF peaks unique to sensitive cells, which suggests that BORIS does not replace CTCF in resistant cells. BORIS preferentially occupied gene regulatory regions-enhancers and promoters (60%)-in resistant cells (Extended Data Fig. 6d, e), which is consistent with its propensity to bind to open chromatin regions²³ (Fig. 2d). Such differential chromatin




Fig. 3 | **BORIS promotes new chromatin interactions in resistant cells. a**, DNA interactions gained in resistant cells based on SMC1A HiChIP analysis. Interaction classes were determined from the genomic locations of the associated anchors (overlapping promoter (Prom) regions (transcription start site (TSS) ± 2 kb), active enhancer (Enh) regions, or CTCF sites only, in that order). Absolute numbers and percentages for each loop type (structural (black), regulatory (blue)) are shown. Cartoon illustrates the spatial proximity induced by DNA looping between these regions. **b**, Fractions of loops bound by BORIS within each interaction class. **c**, Meta-analysis of average CTCF and BORIS ChIP–seq signals in sensitive and resistant cells at the three main anchor types normalized by the number of interactions (n = 2 biological replicates). Anchor sites were centred and extended in both directions (± 2 kb). **d**, ChIP–seq tracks of the indicated proteins in sensitive and resistant cells at the *BORIS* locus

binding at distinct highly expressed genes in resistant versus sensitive cells was commensurate with the MYCN-to-BORIS dependency switch (Extended Data Fig. 6f, g).

The proclivity of aberrantly expressed BORIS for genomic regions associated with active chromatin features in resistant cells suggested that it may, like CTCF and cohesin, regulate gene expression through chromatin looping. Thus, we examined the chromatin looping profiles of sensitive and resistant cells, using cohesin (SMC1A)-based high-throughput chromosome conformation capture followed by chromatin immunoprecipitation (HiChIP)²⁴ (Extended Data Fig. 7a). On the basis of the genomic locations of the associated loop anchors, six classes of interactions were identified²⁵: three longer average interaction loops with a CTCF site on at least one anchor; and three smaller connecting regulatory regions (Fig. 3a, Extended Data Fig. 7b). The overlap of BORIS binding with loop anchors revealed that most (56%) of the 9,487 interactions gained in resistant cells were positive for BORIS (log_2 -transformed fold change > 1; false discovery rate (FDR) < 0.01) (Fig. 3b, Extended Data Fig. 7c). Notably, BORIS was enriched at anchors that were associated with regulatory regions, whereas CTCF binding remained constant, as seen at the BORIS locus itself (Fig. 3c, d). In fact, BORIS binding alone at CTCF-negative loop anchors was sufficient to generate new interactions in resistant cells (Extended Data Fig. 7d).

(representative of two independent experiments), with resistant cell-specific regulatory interactions shown below (HiChIP resistant: pairedend tag (PET) numbers, next to each interaction). Signal intensity is given in the top left corner for each track. **e**, PET interactions in BORIS-depleted (shBORIS) versus control (shCtrl) cells. **f**, Resistant cell-specific loops lost after depletion of BORIS based on loops negative or positive for BORIS binding in shCtrl cells (left), and the odds ratio of losing a loop previously bound by BORIS (right). *P* value determined by two-sided Fisher's exact test. **g**, Meta-analysis of average BORIS, SMC1A and CTCF ChIP-seq signals at resistant cell-specific loop anchors that were lost after depletion of BORIS (n = 2 biological replicates). BORIS depletion at loop anchors inhibits retention of the cohesin complex, and thus prevents the formation of new loops (loop extrusion model). In **a**, **b**, **e** and **f**, n = 3 biological replicates.

To test whether the newly formed interactions in resistant cells were mediated by BORIS binding, we analysed the consequences of BORIS depletion on loop architecture (Extended Data Fig. 7e). Regulatory interactions specific to resistant cells displayed a global shift towards loss after knockdown of BORIS (Fig. 3e), with more than one-quarter of the total interactions lost, of which 63% were positive for BORIS at their anchors (Fig. 3f). Interactions in which anchors were bound by BORIS (especially enhancer–promoter and promoter–promoter interactions) were more likely to be lost after BORIS depletion than those that were not BORIS-bound (Fig. 3f, Extended Data Fig. 7f, g). These results agree with the loop extrusion model²⁶, as BORIS loss resulted in decreased SMC1A binding, preferentially at lost interactions, whereas CTCF binding did not change significantly (Fig. 3g, Extended Data Fig. 7h–j). These data confirm that BORIS is a crucial factor in the looping landscape of resistant cells.

Genes associated with new BORIS-positive regulatory interactions were expressed at higher levels than those associated with BORIS-negative regulatory interactions or genes not associated with new regulatory interactions (Fig. 4a). Because genes that define cell identity are often regulated by super-enhancers in both healthy and cancer cells^{15,27,28}, we characterized the super-enhancer landscape of our cells, observing that the super-enhancers unique to resistant cells were enriched at BORIS-positive regulatory loops (Extended Data Fig. 8a–c).




Fig. 4 | **BORIS-regulated chromatin remodelling supports a phenotypic switch that maintains the resistant state. a**, Left, fold change in expression in counts per million (CPM) of genes involved in resistant cell-specific regulatory interactions that are positive for BORIS binding (n = 1,368) versus those involved in regulatory interactions that are negative for BORIS binding (n = 519) or not associated with a new regulatory interaction (other) (n = 16,151). Centre, fold change in expression of genes involved in resistant cell-specific regulatory interactions positive for BORIS binding and associated with superenhancers (SEs) specific to resistant cells (n = 134) versus those with super-enhancers (n = 720). Right, fold change in expression of genes involved in resistant cell-specific regulatory interactions positive for BORIS binding and associated with resistant cell-specific super-enhancers before and after BORIS knockdown (KD) (n = 134) (P values determined

The presence of such super-enhancers correlated significantly with higher expression of their associated genes in resistant versus sensitive cells (Fig. 4a). These BORIS-positive super-enhancer-associated genes were also enriched for genes that underwent a chromatin state switch from a closed or neutral to an open configuration in resistant cells (Extended Data Fig. 8d, e). Depletion of BORIS resulted in the decreased expression of genes associated with BORIS-positive interactions, especially genes associated with resistant cell-specific superenhancers (Fig. 4a, Extended Data Fig. 8f). These observations suggest that BORIS-mediated alterations in chromatin looping lead to interactions of newly formed super-enhancers with their target genes, which results in their increased expression.

We next sought to identify BORIS-regulated genes that are functionally linked to the resistance phenotype by integrating gene expression, BORIS-mediated looping, super-enhancer landscape and chromatin state. This analysis revealed 89 genes (Supplementary Table), including

by two-sided Wilcoxon rank-sum test). For all box plots, centre lines denote medians; box limits denote twenty-fifth and seventy-fifth percentiles; whiskers denote minima and maxima $(1.5 \times$ the interquartile range). **b**, Highest-ranked transcription factors associated with the resistance phenotype selected based on the presence of at least four of the five indicated features. **c**, ChIP-seq tracks of the indicated proteins in sensitive and resistant cells at the *NEUROG2* locus; regulatory interactions with PET numbers indicated below. **d**, Transcription factor recognition motifs at super-enhancers and promoters (± 2 kb) of the 1,000 highest-expressed genes in resistant and sensitive cells (n = 2 biological replicates) (P values determined by hypergeometric enrichment test). Panels **a**-**c** integrate data of biological replicates from expression microarrays (n = 2), ChIP-seq (n = 2) and HiChIP (n = 3). **e**, Proposed role of BORIS in resistant cells.

13 transcription factors, that are highly expressed during early neural development and are crucial to cell fate decisions^{20,29,30} (Fig. 4b, c, Extended Data Fig. 8g). The expression of these proneural transcription factors paralleled that of BORIS in resistant cells, and was dependent on BORIS-mediated looping, as BORIS depletion led to their downregulation (Extended Data Fig. 8h, i). Moreover, analysis of transcription factor binding sites revealed enrichment of BORIS and several of these proneural transcription factors at the regulatory regions of the highest-expressed genes in resistant cells, whereas sensitive cells were dominated by MYC, MYCN and MAX E-box and E-box-like motifs (Fig. 4d). Similar increased expression of proneural transcription factors with increased BORIS occupancy at their promoters was seen in BORIS-overexpressing E9-resistant SK-N-BE(2) neuroblastoma cells compared with their sensitive counterparts (Extended Data Fig. 8j, k). The high transcriptional activity of these BORIS-regulated genes was also associated with increased binding of the transcriptional activator BRD4, which rendered the resistant cells more sensitive to BET inhibition (Extended Data Fig. 9; Supplementary Note 2). Together, these results indicate the establishment of an alternative transcription factor regulatory network controlled by BORIS-induced chromatin remodelling to support the resistant cell state.

Thus, using a pair of isogenic ALK-inhibitor sensitive and resistant neuroblastoma cell lines, we show that the CTCF paralogue BORIS can promote regulatory DNA interactions that support a phenotypic switch in the context of treatment resistance (Fig. 4e). This mechanism appears relevant to different neuroblastoma cell lines and kinase inhibitors and may extend to other cancers. In Ewing sarcoma, in which overexpression of *BORIS* is associated with metastasis and relapse (Extended Data Fig. 1c), we observed increased BORIS occupancy at regulatory regions in chemotherapy-resistant cell lines (Extended Data Fig. 10; Supplementary Note 3). Further work will establish whether BORIS-mediated alteration of chromatin looping is a general mechanism by which tumour cells co-opt developmental networks to sustain alternative cell states in response to targeted or conventional therapies.

Online content

Any methods, additional references, Nature Research reporting summaries, source data, extended data, supplementary information, acknowledgements, peer review information; details of author contributions and competing interests; and statements of data and code availability are available at https://doi.org/10.1038/s41586-019-1472-0.

Received: 21 July 2017; Accepted: 11 July 2019; Published online: 07 August 2019

- 1. Dixon, J. R. et al. Topological domains in mammalian genomes identified by analysis of chromatin interactions. *Nature* **485**, 376–380 (2012).
- Phillips-Cremins, J. E. et al. Architectural protein subclasses shape 3D organization of genomes during lineage commitment. *Cell* 153, 1281–1295 (2013).
- 3. Hnisz, D. et al. Activation of proto-oncogenes by disruption of chromosome neighborhoods. *Science* **351**, 1454–1458 (2016).
- Katainen, R. et al. CTCF/cohesin-binding sites are frequently mutated in cancer. Nat. Genet. 47, 818–821 (2015).
- Flavahan, W. A. et al. Insulator dysfunction and oncogene activation in *IDH* mutant gliomas. *Nature* 529, 110–114 (2016).
- Loukinov, D. I. et al. BORIS, a novel male germ-line-specific protein associated with epigenetic reprogramming events, shares the same 11-zinc-finger domain with CTCF, the insulator protein involved in reading imprinting marks in the soma. *Proc. Natl Acad. Sci. USA* **99**, 6806–6811 (2002).
- Klenova, E. M., Morse, H. C. III, Ohlsson, R. & Lobanenkov, V. V. The novel BORIS + CTCF gene family is uniquely involved in the epigenetics of normal biology and cancer. Semin. Cancer Biol. 12, 399–414 (2002).
- Martin-Kleiner, I. BORIS in human cancers a review. Eur. J. Cancer 48, 929–935 (2012).
- 9. Garikapati, K. R. et al. Down-regulation of BORIS/CTCFL efficiently regulates cancer stemness and metastasis in *MYCN* amplified neuroblastoma cell line by

modulating Wnt/ β -catenin signaling pathway. Biochem. Biophys. Res. Commun. **484**, 93–99 (2017).

- Cheung, N. K. & Dyer, M. A. Neuroblastoma: developmental biology, cancer genomics and immunotherapy. *Nat. Rev. Cancer* 13, 397–411 (2013).
- Monk, M., Hitchins, M. & Hawes, S. Differential expression of the embryo/cancer gene ECSA(DPPA2), the cancer/testis gene BORIS and the pluripotency structural gene OCT4, in human preimplantation development. *Mol. Hum. Reprod.* 14, 347–355 (2008).
- Brodeur, G. M., Seeger, R. C., Schwab, M., Varmus, H. E. & Bishop, J. M. Amplification of N-myc in untreated human neuroblastomas correlates with advanced disease stage. *Science* 224, 1121–1124 (1984).
- George, R. E. et al. Activating mutations in ALK provide a therapeutic target in neuroblastoma. *Nature* 455, 975–978 (2008).
- Galkin, A. V. et al. Identification of NVP-TAE684, a potent, selective, and efficacious inhibitor of NPM-ALK. *Proc. Natl Acad. Sci. USA* 104, 270–275 (2007).
- Chipumuro, E. et al. CDK7 inhibition suppresses super-enhancer-linked oncogenic transcription in MYCN-driven cancer. Cell 159, 1126–1139 (2014).
- Durbin, A. D. et al. Selective gene dependencies in MYCN-amplified neuroblastoma include the core transcriptional regulatory circuitry. Nat. Genet. 50, 1240–1246 (2018).
- van der Maaten, L. & Hinton, G. Visualizing data using t-SNE. J. Mach. Learn. Res. 9, 2579–2605 (2008).
- Zeid, R. et al. Enhancer invasion shapes MYCN-dependent transcriptional amplification in neuroblastoma. *Nat. Genet.* 50, 515–523 (2018).
- Muhar, M. et al. SLAM-seq defines direct gene-regulatory functions of the BRD4-MYC axis. Science 360, 800–805 (2018).
- Wegner, M. & Stolt, C. C. From stem cells to neurons and glia: a Soxist's view of neural development. *Trends Neurosci.* 28, 583–588 (2005).
- Gao, Y. et al. Overcoming resistance to the THZ series of covalent transcriptional CDK inhibitors. *Cell Chem. Biol.* 25, 135–142 (2018).
- Pugacheva, E. M. et al. Comparative analyses of CTCF and BORIS occupancies uncover two distinct classes of CTCF binding genomic regions. *Genome Biol.* 16, 161 (2015).
- 23. Bergmaier, P. et al. Choice of binding sites for CTCFL compared to CTCF is driven by chromatin and by sequence preference. *Nucleic Acids Res.* **46**, 7097–7107 (2018).
- Mumbach, M. R. et al. HiChIP: efficient and sensitive analysis of protein-directed genome architecture. Nat. Methods 13, 919–922 (2016).
- Dowen, J. M. et al. Control of cell identity genes occurs in insulated neighborhoods in mammalian chromosomes. *Cell* 159, 374–387 (2014).
- Sanborn, A. L. et al. Chromatin extrusion explains key features of loop and domain formation in wild-type and engineered genomes. *Proc. Natl Acad. Sci.* USA 112, E6456–E6465 (2015).
- Whyte, W. A. et al. Master transcription factors and mediator establish super-enhancers at key cell identity genes. *Cell* 153, 307–319 (2013).
- Hnisz, D. et al. Super-enhancers in the control of cell identity and disease. Cell 155, 934–947 (2013).
- Kumar, J. P. The sine oculis homeobox (SIX) family of transcription factors as regulators of development and disease. *Cell. Mol. Life Sci.* 66, 565–583 (2009).
- Dennis, D. J., Han, S. & Schuurmans, C. bHLH transcription factors in neural development, disease, and reprogramming. *Brain Res.* 1705, 48–65 (2019).

Publisher's note: Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

© The Author(s), under exclusive licence to Springer Nature Limited 2019

METHODS

Cell lines. Human neuroblastoma cell lines Kelly and SK-N-BE(2) and human Ewing sarcoma cell lines TC-32, TC-71 and CHLA-10^{31,32} were obtained from the Children's Oncology Group cell line bank. Human neuroblastoma cell line SK-N-SH and human embryonic kidney cell line HEK293T were obtained from the American Type Culture Collection. Cell line authenticity was confirmed by genotyping, and cells were tested negative for mycoplasma contamination every 3 months. All cells except HEK293T were grown in RPMI-1640 medium supplemented with 10% fetal bovine serum (FBS) and 1% penicillin/streptomycin (Life Technologies). HEK293T cells were grown in DMEM medium supplemented with 10% FBS and 1% penicillin/streptomycin (Life Technologies). Resistant cells were grown in the presence of either the ALK inhibitor, TAE684¹⁴ (Kelly and SK-N-SH) or the CDK12 inhibitor, E9²¹ (SK-N-BE(2)).

Compounds. TAE684 and E9 were synthesized in-house in the Gray laboratory and JQ1³³ was obtained from J.Qi's laboratory at the Dana-Farber Cancer Institute (DFCI). Ceritinib³⁴, lorlatinib³⁵ and I-BET726³⁶ were purchased from Selleck Chemicals.

Synthetic RNA spike-in and microarray analysis. Total RNA and sample preparation was performed as previously described³⁷. In brief, cells were either incubated in medium containing DMSO, TAE684 (1 μ M) or JQ1 (2.5 μ M), or infected with shRNA (Ctrl or BORIS) for 24 h. Cell numbers were determined using a Countess II cell counter (Life Technologies) before lysis and RNA extraction. Biological duplicates (equivalent to 5 × 10⁶ cells per replicate) were collected and homogenized in 1 ml of TRIzol Reagent (Ambion), purified using the mirVANA miRNA isolation kit (Ambion) following the manufacturer's instructions and re-suspended in 50 μ nuclease-free water (Ambion). Total RNA was spiked-in with ERCC RNA Spike-In Mix (Ambion), treated with DNA-free DNase I (Ambion) and analysed on an Agilent 2100 Bioanalyzer (Agilent Technologies) for integrity. RNA with the RNA Integrity Number above 9.8 was hybridized to Affymetrix GeneChip PrimeView Human Gene Expression arrays (Affymetrix).

Antibodies. The following antibodies were used: N-MYC (9405), N-MYC (51705), cleaved PARP (9541), cleaved caspase 3 (9661), ALK (3333), AKT (4691), pAKT-T308 (9275), pAKT-S473 (#9271), ERK (4695), pERK (4377), S6 (2217), pS6 (4857), STAT3 (4904), pSTAT3 (9131), ABCB1 (12683), SOX2 (3579), β -actin (4967), CTCF (3417), normal rabbit IgG (2729) and HRP anti-mouse IgG (7076) from Cell Signaling Technology; HRP anti-rabbit IgG (sc-2357) from Santa Cruz Biotechnology; BRD4 (A301-985A100) and SMC1A (A300-055A) from Bethyl Laboratories; CTCF (07-729), SOX9 (AB5535) and H3K27me3 (07-449) from Millipore; pALK-Y1507 (ab73996), BORIS (ab187163) and H3K27ac (ab2729) from Abcam; BORIS (NBP2-52405) from NOVUS Biologicals; BORIS (39851) from Active Motif; SIX1 (HPA001893) from Sigma-Aldrich; and Vysis LSI N-MYC (2p24) SpectrumGreen/Vysis CEP 2 SpectrumOrange Probe (07J72-001) from Abbott.

Cell viability and growth curve assays. Viability and growth experiments were performed using the CellTiter-Glo Luminescent Cell Viability Assay (Promega) according to the manufacturer's instructions, as previously described³⁸. Cells were plated in 96-well plates at a seeding density of 4×10^3 cells per well. For growth assays, the cells were analysed each day until day 5. For viability, after 24 h, the cells were treated with various concentrations of the indicated drug (ranging from 1 nM to 10 μ M except for I-BET726: 2 nM to 20 μ M). DMSO without drug served as a negative control. After 72 h of incubation, cells were analysed for cell viability and IC₅₀ values were determined using a nonlinear regression curve fit with GraphPad Prism 6 software.

Cell-cycle analysis. Cell-cycle analysis was performed 24 h after cell plating using propidium iodide staining, as previously described¹⁵. Cells fixed with 80% ethanol overnight at 4°C were resuspended in PBS supplemented with 0.1% Triton X-100 (Sigma-Aldrich), 25 mg ml⁻¹ propidium iodide (BD Biosciences) and 0.2 mg ml⁻¹ RNase A (Sigma-Aldrich). After 45 min at 37 °C in the dark, analysis was performed on a FACSCalibur flow cytometer (BD Biosciences). Cell-cycle profiles were plotted as histograms generated using FlowJo software (FLOWJO).

Western blotting. Cell or tumour tissue was lysed in NP-40 buffer (Invitrogen) containing a 1× complete protease inhibitor tablet (Roche) per 10 ml buffer and a cocktail of phosphatase inhibitors (Roche). Protein concentration was measured using the DC Protein Assay (Bio-Rad); protein (50 μ g) was denatured in LDS sample buffer with reducing agent (Invitrogen), separated on precast 4–12% Bis-Tris gels (Life Technologies) and transferred to nitrocellulose membranes (Bio-Rad). Membranes were incubated in blocking buffer (5% dry milk in TBS with 0.2% Tween-20) for 1 h, and then incubated in the primary antibody in blocking buffer overnight at 4°C. Chemiluminescent detection was performed with the appropriate secondary antibodies and developed using Genemate Blue ultra-autoradiography film (VWR). The actin loading controls for the protein samples shown in the immunoblots of the following panels (two independent experiments) are the same because the samples were run on a single gel but probed for pALK, ALK (Extended Data

Fig. 2a), MYCN (Extended Data Fig. 3e) and BORIS (Extended Data Fig. 4a), respectively.

Co-immunoprecipitation. Cells were collected in immunoprecipitation lysis buffer (50 mM Tris-HCl buffer (pH 7.4), 100 mM NaCl, 1% Triton-100, 1 mM PMSF), containing a 1× complete protease inhibitor tablet (Roche) per 10 ml buffer and a cocktail of phosphatase inhibitors (Roche). Homogenates were centrifuged at 20,000g for 10 min at 4°C to obtain supernatants. DNase I (approximately 1 U ml⁻¹) was used to degrade DNA in supernatants by incubation for 1 h at room temperature. Co-immunoprecipitation of endogenously expressed proteins was performed using protein A Dynabeads (Invitrogen), according to the manufacturer's instructions. In brief, antibody-conjugated Dynabeads were incubated with purified cell lysates to immunoprecipitate the target antigen. Antibodies used for immunoprecipitation were CTCF (3417, Cell Signaling Technology) and BORIS (NBP2-52405, NOVUS Biologicals). The elution step was conducted by heating the beads for 10 min at 95 °C in lithium dodecyl sulfate (LDS) sample buffer with reducing agent (Invitrogen), after which western blotting was performed using the following antibodies: CTCF (3417, Cell Signaling Technology) and BORIS (9851, Active Motif).

Plasmids, shRNA knockdown and overexpression systems. pLKO.1 shRNA constructs (control: SHC007; MYCN: 1-TRCN0000020694 and 2-TRCN0000363425; BORIS: 3-TRCN0000370229 and 4-TRCN0000365141; BRD4: A-TRCN0000318771 and B-TRCN0000196576) were purchased from Sigma-Aldrich and pLKO.1 GFP shRNA was a gift from D. Sabatini (Addgene plasmid 30323)³⁹. Overexpression constructs were generated by cloning BORIS cDNA into the Tet-inducible pInducer20 vector, provided by S. Elledge (Addgene plasmid 44012)⁴⁰. Production of lentiviral particles and subsequent infection were performed as previously described³⁸. The lentivirus was packaged by co-transfection of either pLKO.1 or pInducer20 plasmid with the helper plasmids, pCMV-deltaR8.91 and pMD2.G-VSV-G into HEK293T cells using TransIT-LT1 Transfection Reagent (Mirus Bio LLC). Virus-containing supernatants were collected 48 h after transfection. Cells were infected with 8 µg ml⁻¹ polybrene (Sigma-Aldrich) and 24-48 h later selected with puromycin (pLKO.1) (Sigma-Aldrich) and then collected at appropriate time points. When using the Tet-inducible system for BORIS overexpression, induction of gene expression was achieved by treating cells every 2–3 days with doxycycline $(0.2 \ \mu g \ ml^{-1})$ for a total duration of 37 days.

qRT-PCR. RNA isolation and PCR amplification were performed as previously described³⁸, except that the RT-PCR was performed using the SuperScript III First-Strand system (Life Technologies). Total RNA was isolated from cell lines with the RNeasy kit (Qiagen). One microgram of purified RNA was reverse transcribed using Superscript III First-Strand (Invitrogen) according to the manufacturer's protocol, and quantitative PCR was performed using SYBR Green on a Viia7 Real-Time PCR system (Thermo Fisher Scientific). All experiments were performed in biological triplicates unless stated otherwise. Each individual biological sample was amplified by qPCR in technical replicates and normalized to actin as an internal control. Amplification was carried out with primers specific to the genes to be quantified (sequences available on request).

Sequence analysis. The kinase domain of ALK was amplified from cDNA extracted from sensitive and resistant cells using the HotStar HiFidelity Polymerase Kit (Qiagen). The PCR products were cloned into the pGEM-T vector (Promega) and confirmed by sequencing.

RTK array. The Human Phospho-RTK Array Kit (R&D Systems) was used as previously described³⁸. Cell lysate (500 μ g) was incubated on a phospho-RTK membrane array (ARY001B) according to the manufacturer's instructions. Target proteins were captured with their respective antibodies. After washing, the proteins were incubated with a phosphotyrosine antibody conjugated to horseradish peroxidase to allow the detection of captured phosphorylated RTKs.

Fluorescent in situ hybridization. Fluorescent in situ hybridization (FISH) analyses were performed using a Vysis LSI N-MYC (2p24) SpectrumGreen/Vysis CEP 2 SpectrumOrange Probe (Vysis), in accordance with the manufacturer's instructions.

Immunohistochemistry. All human tumour specimens (formalin-fixed paraffin-embedded slides) were obtained under an Institutional Review Board (IRB)-approved protocol of the Dana-Farber/Boston Children's Cancer and Blood Disorders Center, and informed consent was obtained from all subjects. Staining was performed by Applied Pathology Systems using the ImmPRESS Excel Amplified HRP Polymer Staining Kit (MP-7601, Vector Laboratories) on a Dako Autostainer (Agilent Technologies). Sections were deparaffinized, rehydrated, and subjected to antigen retrieval in citrate-based buffer on a steamer for 25 min. Slides were blocked with BLOXALL blocking solution and 2.5% horse serum sequentially before a 1-h incubation with BORIS antibody at 1:50 dilution (ab187163, Abcam). Sections were then incubated with anti-rabbit amplifier antibody and ImmPRESS Excel Amplified HRP Polymer Reagent sequentially before incubation with ImmPACT DAB EqV Substrate. Finally, slides were counterstained with haematoxylin, followed by dehydration and the addition of coverslips. **Bisulfite sequencing.** Methylation analysis of *BORIS* (NCBI RefSeq NC_000020.11, spanning nucleotides chr20: 57,524,203–57,525,234 on GRCh38.p7 assembly) was performed using a bisulfite sequencing assay. Genomic DNA (500 ng) was treated with the EZ DNA Methylation-Lightning Kit (Zymo Research), followed by PCR using ZymoTaq Polymerase premix (Zymo Research) and specific primers designed using the Zymo bisulfite primer seeker (http://www.zymore-search.com/tools/bisulfite-primer-seeker/; sequences available on request). PCR products were then sequenced for the assessment of CpG site-specific DNA methylation in the *BORIS* promoter region.

Growth assay. After shRNA-mediated knockdown of BORIS, cells were reseeded at a density of 4×10^5 cells per well in 6-well plates. At 48 and 120 h of incubation, cells were stained with trypan blue (Sigma-Aldrich) and counted on a Countess II cell counter (Life Technologies).

Mouse experiments. All mouse experiments were performed with approval from the Institutional Animal Care and Use Committee (IACUC) of the DFCI. Three mouse experiments were performed: (i) to assess the tumorigenic potential of resistant cells in vivo; (ii) to assess that resistance to TAE684 was maintained in vivo; and (iii) to assess the effect of JQ1 on resistant cells in vivo. All experiments were performed using subcutaneous cell xenograft models generated by injecting 2×10^6 sensitive or resistant Kelly neuroblastoma cells into the flanks of NU/NU (Crl:NU-Foxn1nu) (Charles River Laboratories) or NU/NU (CrTac:NCr-Foxn1nu) (Taconic) 7-week-old female mice. Mice were randomized into groups of equal average volumes, and investigators were not blinded to group allocation during data collection. (i) To assess the tumorigenic potential of resistant cells in absence of treatment, mice with established disease (mean tumour volume of 200 mm³) were monitored for up to 23 days (n = 4 per group). Tumours were obtained, dissociated and used to establish cell lines and for assessment of mRNA levels, protein expression and sensitivity to TAE684. (ii) To ensure that the in vitro resistance to TAE684 was maintained in vivo, mice with established disease were divided into two cohorts and were treated with either TAE684 (10 mg kg⁻¹) or vehicle control by oral gavage once daily (n = 8per group), and were monitored for up to 56 days from start of treatment. (iii) To assess the sensitivity of resistant cells to BRD4 inhibition, mice with established disease were divided into two cohorts and treated with either JQ1 (50 mg kg⁻¹) or vehicle control intraperitoneally (i.p.) once daily (n = 6per group), and were monitored for up to 87 days from start of treatment. For all experiments, disease burden was quantified using electronic caliper measurements (2-3 times a week) and mouse weights were monitored at least twice a week. Tumour volumes were calculated using the modified ellipsoid formula⁴¹: $\frac{1}{2}(\text{length} \times \text{width}^2)$. Animals were euthanized when tumour volumes reached 1,500-2,000 mm³ based on institutional IACUC criteria for maximum tumour volumes. In none of the experiments were the institutional limits for tumour volumes (<2,000 mm³ measurement preceding the day of euthanization) exceeded.

ChIP-seq. ChIP was carried out as previously described¹⁵ with minor changes as described. Approximately 1×10^7 cells were crosslinked for 10 min at room temperature with 1% formaldehyde (Thermo Scientific) in PBS followed by quenching with 0.125 M glycine for 5 min. The cells were then washed twice in ice-cold PBS, and the cell pellets flash frozen and stored at -80 °C. Fifty microlitres of protein G Dynabeads per sample (Invitrogen) were blocked with 0.02% Tween20 (w/v) in PBS. Magnetic beads were loaded with 10 μg each of antibody and incubated overnight at 4 °C. Crosslinked cells were lysed, placed in sonication buffer with 0.2% SDS, placed on ice and chromatin was sheared using a Misonix 3000 sonicator (Misonix) at the following settings: 10 cycles, each for 30 s on, followed by 1 min off, at a power of approximately 20 W. The lysates were then centrifuged for 10 min at 4°C, supernatants collected and diluted with an equal amount of sonication buffer to reach a final concentration of 0.1% SDS. The sonicated lysates were incubated overnight at 4°C with the antibody-bound magnetic beads, washed with low-salt buffer (50 mM HEPES-KOH (pH 7.5), 0.1% SDS, 1% Triton X-100, 0.1% sodium deoxycholate, 1 mM EGTA, 1 mM EDTA, 140 mM NaCl and 1 \times complete protease inhibitor), high-salt buffer (50 mM HEPES-KOH (pH 7.5), 0.1% SDS, 1% Triton X-100, 0.1% sodium deoxycholate, 1 mM EGTA, 1 mM EDTA, 500 mM NaCl and 1× complete protease inhibitor), LiCl buffer (20 mM Tris-HCl (pH 8), 0.5% NP-40, 0.5% sodium deoxycholate, 1 mM EDTA, 250 mM LiCl and $1 \times$ complete protease inhibitor) and Tris-EDTA buffer. DNA was then eluted in elution buffer (50 mM Tris-HCl (pH 8.0), 10 mM EDTA, 1% SDS), and highspeed centrifugation was performed to pellet the magnetic beads and collect the supernatants. The crosslinking was reversed overnight at 65 °C. RNA and protein were digested using RNase A and proteinase K, respectively, and DNA was purified with phenol chloroform extraction and ethanol precipitation. Purified ChIP DNA was used to prepare Illumina multiplexed sequencing libraries using the NEBNext Ultra II DNA Library Prep kit and the NEBNext Multiplex Oligos for Illumina

(New England Biolabs) according to the manufacturer's protocol. Libraries with distinct indexes were multiplexed and run together on the Illumina NextSeq 500 (SY-415-1001, Illumina) for 75 bases in single-read mode.

HiChIP. HiChIP was performed as previously described²⁴ with a few modifications. Approximately 1×10^7 cells were crosslinked for 10 min at room temperature with 1% formaldehyde in growth medium and quenched in 0.125 M glycine. After washing twice with ice-cold PBS, the supernatant was aspirated and the cell pellet flash frozen in liquid nitrogen. Crosslinked cell pellets were thawed on ice, resuspended in 1 ml of ice-cold Hi-C lysis buffer (10 mM Tris-HCl (pH 8.0), 10 mM NaCl, 0.2% NP-40 and $1 \times$ complete protease inhibitor) and incubated at 4°C for 30 min with rotation. Nuclei were pelleted by centrifugation for 5 min at 4°C and washed once with 500 µl of ice-cold Hi-C lysis buffer. After removing the supernatant, nuclei were resuspended in 100 μl of 0.5% SDS and incubated at 62 °C for 10 min. SDS was quenched by adding $335\,\mu l$ of 1.5% Triton X-100 and incubating for 15 min at 37 °C. After the addition of 50 μ l of 10 \times NEB Buffer 2 (New England Biolabs, B7002) and 375 U of MboI restriction enzyme (New England Biolabs, R0147), chromatin was digested at 37 °C for 2 h with rotation. After digestion, MboI enzyme was heat-inactivated by incubating the nuclei at 62 °C for 20 min. To fill in the restriction fragment overhangs and mark the DNA ends with biotin, 52 µl of fill-in master mix, containing 37.5 µl of 0.4 mM biotin-dATP (Invitrogen, 19524016), 1.5 µl of 10 mM dCTP (Invitrogen, 18253013), 1.5 µl of 10 mM dGTP (Invitrogen, 18254011), 1.5 µl of 10 mM dTTP (Invitrogen, 18255018), and 10 μ l of 5 U μ l⁻¹ DNA Polymerase I, Large (Klenow) Fragment (New England Biolabs, M0210), were added and the tubes were incubated at 37 °C for 1 h with rotation. Proximity ligation was performed by the addition of 948 μ l of ligation master mix, containing 150 μ l of 10 \times NEB T4 DNA ligase buffer (New England Biolabs, B0202), 125 μ l of 10% Triton X-100, 7.5 μ l of 20 mg ml⁻¹ BSA (New England Biolabs, B9000), 10 μl of 400 U μl^{-1} T4 DNA ligase (New England Biolabs, M0202), and 655.5 μ l of water, and incubation at room temperature for 4 h with rotation. After proximity ligation, nuclei were pelleted by centrifugation for 5 min and resuspended in 1 ml of ChIP sonication buffer (50 mM HEPES-KOH (pH 7.5), 140 mM NaCl, 1 mM EDTA pH 8.0, 1 mM EGTA (pH 8.0), 1% Triton X-100, 0.1% sodium deoxycholate, 0.1% SDS and $1 \times$ complete protease inhibitor). Nuclei were sonicated using a Misonix 3000 sonicator (Misonix) at the following settings: 12 cycles, each for 30 s on, followed by 1 min off, at a power of approximately 20 W. Sonicated chromatin was clarified by centrifugation for 15 min at 4°C and the supernatant was transferred to a tube. Sixty microlitres of protein G Dynabeads (Invitrogen) were washed three times and resuspended in 50 µl sonication buffer. Washed beads were then added to the sonicated chromatin and incubated for 1 h at 4 °C with rotation. Beads were then separated on a magnetic stand and the supernatant was transferred to a new tube. Seventy-five microlitres of protein G Dynabeads pre-incubated overnight at 4 °C with 10 µg of anti-SMC1A antibody (Bethyl A300-055A) or 10 µg of BORIS antibody (Abcam, ab187163) were added to the tube and incubated overnight at 4°C with rotation. Beads were then separated on a magnetic stand and washed twice with 1 ml of sonication buffer, followed by once with 1 ml high-salt sonication buffer (50 mM HEPES-KOH (pH 7.5), 500 mM NaCl, 1 mM EDTA pH 8.0, 1 mM EGTA (pH 8.0), 1% Triton X-100, 0.1% sodium deoxycholate, 0.1% SDS), once with 1 ml of LiCl wash buffer (20 mM Tris-HCl (pH 8.0), 1 mM EDTA pH 8.0, 250 mM LiCl, 0.5% NP-40, 0.5% sodium deoxycholate, 0.1% SDS) and once with 1 ml of TE buffer with salt (10 mM Tris-HCl (pH 8.0), 1 mM EDTA pH 8.0, 50 mM NaCl). Beads were then resuspended in 200 μ l of elution buffer (50 mM Tris-HCl (pH 8.0), 10 mM EDTA pH 8.0, 1% SDS) and incubated at 65 °C for 15 min. To purify the eluted DNA, RNA was degraded by the addition of 8.5 µl of 10 mg ml⁻¹ RNase A and incubation at 37 °C for 2 h. Protein was degraded by the addition of 20 μ l of 10 mg ml⁻¹ proteinase K and incubation at 55 °C for 45 min. Samples were then incubated at 65 °C overnight to reverse crosslink protein-DNA complexes. DNA was then purified using Zymo ChIP DNA Clean and Concentrator columns (Zymo, D5205) according to the manufacturer's protocol and eluted in 14 µl water. The amount of eluted DNA was quantified by Qubit dsDNA HS kit (Invitrogen, Q32854). Tagmentation of ChIP DNA was performed using the Illumina Nextera DNA Library Prep Kit (Illumina, FC-121-1030). First, 5 µl of MyOne Streptavidin C1 Dynabeads (Invitrogen, 65001) was washed with 1 ml of Tween wash buffer (5 mM Tris-HCl (pH 7.5), 0.5 mM EDTA (pH 8.0), 1 M NaCl, 0.05% Tween-20) and resuspended in 10 μ l of 2× biotin binding buffer (10 mM Tris-HCl (pH 7.5), 1 mM EDTA pH 8.0, 2 M NaCl). Then, 25 ng of purified DNA was added in a total volume of 10 μ l water to the beads and incubated at room temperature for 15 min with agitation every 5 min. After capture, beads were separated with a magnet and the supernatant was discarded. Beads were then washed twice with 500 µl of Tween wash buffer, incubating at 55 °C for 2 min with shaking for each wash. Beads were resuspended in 25 μl of Nextera Tagment DNA buffer. To tagment the captured DNA, 1 µl of Nextera Tagment DNA Enzyme 1 was added with 24 µl of Nextera Resuspension Buffer and samples were incubated at 55 °C for 10 min with shaking. Beads were separated on a magnet and supernatant was discarded. Beads were

washed twice with 500 µl of 50 mM EDTA at 50 °C for 30 min, washed twice with 500 µl of Tween wash buffer at 55 °C for 2 min each, and finally washed once with 500 µl of 10 mM Tris-HCl (pH 7.5) for 1 min at room temperature. Beads were separated on a magnet and supernatant was discarded. To generate the sequencing library, PCR amplification of the tagmented DNA was performed while the DNA was still bound to the beads. Beads were resuspended in 15 µl of Nextera PCR Master Mix, 5 µl of Nextera PCR Primer Cocktail, 5 µl of Nextera Index Primer 1, 5 µl of Nextera Index Primer 2 and 20 µl water. DNA was amplified with 9–10 cycles of PCR. After PCR, beads were separated on a magnet and the supernatant containing the PCR-amplified library was transferred to a new tube, purified using Zymo DNA Clean and Concentrator columns (Zymo, D5205) according to the manufacturer's protocol and subjected to 2 \times 100 paired-end sequencing to the manufacturer's protocol and subjected to 2 \times 100 paired-end sequencing using an Illumina HiSeq 2500 system (SY–401–2501, Illumina).

scRNA-seq. Kelly cells (sensitive, intermediate and resistant states) were grown to 70% confluence in T75 culture flasks. In brief, growth medium was aspirated and cells were treated with 0.25% Trypsin/EDTA for 3 min at 37 °C, after which cells were washed twice with $1 \times PBS$. Cells were then resuspended into single cells at a concentration of 1×10^6 per ml in $1 \times$ PBS with 0.4% BSA for 10x Genomics processing. The sorted cell suspensions were loaded onto a 10x Genomics Chromium instrument to generate single-cell gel beads in emulsion (GEMs). Approximately 5,000 cells were loaded per channel. scRNA-seq libraries were prepared using the following Single Cell 3' Reagent Kits: Chromium Single Cell 3' Library & Gel Bead Kit v2 (PN-120237), Single Cell 3' Chip Kit v2 (PN-120236) and i7 Multiplex Kit (PN-120262) (10x Genomics) as previously described⁴², and following the Single Cell 3' Reagent Kits v2 User Guide (Manual Part CG00052 Rev A). Libraries were run on an Illumina HiSeq 4000 system (SY-401-4001, Illumina) as 2×150 pairedend reads, one full lane per sample, for approximately >90% sequencing saturation. Genomics analysis: direct comparison of CTCF and BORIS expression in healthy and tumour samples. To assess the expression levels and range of BORIS and CTCF in healthy and tumour cells all GTEx, TCGA and TARGET datasets were downloaded and converted to FPKM values and displayed as [log2(FPKM + 1)] (Extended Data Fig. 1a, b).

Association of BORIS with prognostic features. For each dataset, processed values were extracted from the Gene Expression Omnibus (GEO) and scaled values were created by normalizing the expression levels by the minimum mean value of the conditions that were compared, Esi,j = Ei,j/min(average(Ej)). The two-sided Wilcoxon rank-sum test on the original values was used to determine statistical differences between the compared conditions (Extended Data Fig. 1c and Extended Data Fig. 4f).

Microarray data analysis. Microarray data were analysed using a custom CDF file (GPL16043) that contained the mapping information of the ERCC probes used in the spike-in RNAs. The arrays were normalized as previously described³⁷. In brief, all microarray chip data were imported in R (https://www.r-project.org/, v.3.1.3) using the affy package⁴³ (v.1.44.0), converted into expression values using the expresso command, normalized to take into account the different numbers of cells and spike-ins used in the different experiments and renormalized using loess regression fitted to the spike-in probes. Sets of differentially expressed genes were obtained using the limma package⁴⁴ (v.3.22.7) and a FDR value of 0.05. Spike-in normalized absolute expression values (counts) were normalized to CPM as a measurement of relative gene expression concentrations per condition. Total number of transcripts per sample was determined as the total number of counts after spike-in normalization and the *BORIS* shRNA sample was first normalized to the transfection protocol.

ChIP-seq analysis. For all ChIP-seq samples, high-quality data were confirmed using the Fastqc tool (v.0.11.5) and samples were aligned to the human genome (build hg19, GRCh37.75) with STAR (v.2.5.1b_modified) and the parameters 'alignIntronMax 1-alignEndsType EndToEnd-outFilterMultimapNmax 1-out-FilterMismatchMax 5'. Next, non-duplicate reads that mapped to the reference chromosomes were retained using Samtools (v.1.3.1) and MarkDuplicates (v.2.1.1) from Picard tools. For each experimental replicate, antibody enrichment was assessed using the plotFingerprint command from deepTools (v.2.2.4). Peaks were identified with MACS2 (v.2.1.1) for narrow peaks (BORIS, CTCF, BRD4, Pol2, MYCN) with the parameters '-q 0.01-call-summits' and for broad peaks (H3K27ac, H3K27me3) with the parameters '-broad-cutoff 0.01'. Peaks overlapping regions with known artefact regions (http://mitra.stanford.edu/kundaje/ akundaje/release/blacklists/) were blacklisted out. Input normalized bedgraph tracks were created with the deepTools command bamCompare and the param $eters\ `-scaleFactorsMethod=readCount-ratio=subtract-binSize=50-number-inter$ OfProcessors=4-extendReads=200'. Subsequently, negative values were set to zero and counts were scaled to RPM per bp to account for differences in library size. Bigwig files were created with bedGraphToBigWig (v.4). ChIP-seq replicates (n = 2) were merged at the BAM level after assessment of strong correlation with the deepTools command 'multiBigwigSummary BED-file' using all replicate bigwigs and identified peaks as input. Identification of peaks and generation of tracks were then repeated for these merged files and used for further analyses. Downstream analyses for ChIP-seq and other genomic interval data was performed in R (https://www.r-project.org/) (v.3.5.1) using the data.table (v.1.12.2) package.

Gencode annotation and isoform selection. Gencode (http://www.gencodegenes. org/, release 19) annotation was used and for each gene the most likely isoform was selected based on data-driven criteria. In brief, only genes that were part of the Refseq transcriptome annotation and with a minimum length of 1 kb were considered. Next, isoforms were prioritized according to increased deposition of Pol2 and H3K27ac reads on the TSS, transcript length and alphabet rank, in that order, until only one transcript was selected for each gene.

Cell-type-specific binding patterns. To determine the cell specificities of BORIS and CTCF peaks, we first combined all peaks identified by MACS2 and merged the peak regions that overlapped by at least 50%. A 50% threshold was empirically selected to avoid merging peaks that had clear and distinct summits. Next, normalized BORIS or CTCF read densities were calculated for each region and a ratio [log2(resistant/sensitive)] was calculated. Peak regions with a twofold density increase or decrease were classified as resistant or sensitive cell-specific peaks, respectively, whereas other regions were denoted as 'shared' to indicate that these peaks had similar BORIS or CTCF deposition in both cell types (Fig. 2a, b and Extended Data Fig. 6a). To explore the proximity of BORIS and CTCF peaks and how they were altered during the transition from sensitive to resistant cells, we overlapped all shared and cell-type-specific peaks from both cell types in the least stringent way (minimum 1-bp overlap) (Fig. 2c and Extended Data Fig. 6a).

Genomic enrichment of peak-binding sites. To identify genomic locations with BORIS or CTCF binding we determined the number of peaks that overlapped with at least 25% of known functional regions in the following order: (i) broad promoter (±2 kb TSS); (ii) BRD4⁺ H3K27ac⁺ (active) enhancers; (iii) BRD4⁻ H3K27ac⁺ enhancers; (iv) exons; (v) introns; (vi) repressed chromatin represented by H3K27me3 broad peaks; or (vii) other (if the peak was outside the aforementioned regions) (Extended Data Fig. 6d). Enrichment of ChIP-seq binding at resistant cell BORIS peaks was performed by extending BORIS summits by 1 kb in both directions and calculating the normalized read densities in 50-bp bins (Fig. 2d). Genomic enrichment of regulatory regions. To visualize the enrichment of CTCF and BORIS at regulatory regions (enhancers and promoters) and the differences between sensitive and resistant cells, a metagene analysis for CTCF and BORIS occupancies was performed for all H3K27ac enhancer regions and gene promoters. All TSSs were extended in both directions by 2 kb and binned in 50-bp bins, and each enhancer (start-end) was divided into 40 equally sized bins and extended with 2 kb in both directions and these extended regions were binned in 50-bp bins. Normalized bedgraph files were used to calculate read density (RPM per bp). An aggregated summary profile was created for each cell type. To account for different numbers of identified enhancers in both cells types we calculated a normalization factor (N resistant enhancers/N sensitive enhancers) to divide each aggregated read density (Extended Data Fig. 6e).

HiChIP processing and quality control. For all SMC1A-based HiChIP datasets, raw reads were first trimmed to a uniform length of 50 bp using trimmomatic⁴⁵(v.0.36) and were then processed using the HiC-Pro (v.2.10.0) pipeline⁴⁶ with default settings for the human genome (build hg19, GRCh37.75) and corresponding MboI cut sites. To perform intra- and inter-correlation analysis for biological replicates, forward and reverse reads from the HiC-Pro output were merged together to generate one-dimensional SMC1A BAM profiles. Genome-wide Spearman correlation in 5-kb bins was computed for all merged genomic anchor regions on those merged BAMs for all replicates using the 'multiBamSummary BED-file' command from deepTools (Extended Data Fig. 7a, e).

HiChIP loop calling and differential looping analysis. Loops were directly called from the HiC-Pro output using hichipper⁴⁷ (v.0.7.3), with parameter 'peaks = combined, all', and subsequently diffloop⁴⁷ (v.1.10.0) with default settings. Only loops that were detected in all three biological replicates of a sample (sensitive, resistant, shCtrl or shBORIS) with a minimum of five paired-end tags in total and an FDR < 0.01 were retained for further analysis. To call differential loops between samples, the quickAssocVoom function was used and significantly different loops were either considered reinforced (mango.FDR < 0.01 and log₂-transformed fold change > 1) or lost (mango.FDR < 0.01 and log₂-transformed fold change < -1). Classification of HiChIP interactions. SMC1A-based HiChIP interactions (loops) were classified as previously described⁴⁸ with minor adaptations. Associated anchors of loops were overlapped with our ChIP-seq peaks (CTCF, BORIS, H3K27ac, BRD4) and promoter regions (TSS \pm 2 kb), requiring a minimum 1-bp overlap. Each anchor was then independently classified according to its overlap profile, following a hierarchical tree. If an anchor overlapped a promoter, an enhancer (BRD4 + H3K27ac), or a CTCF peak, it was classified as

promoter-, enhancer- or CTCF-anchor, in that order. If there was no overlap, the anchor was considered 'other'. By combining these four anchor classes we discriminated 10 different interaction classes. We excluded from further analyses any interaction that contained an anchor classified as other, which also represented on average much shorter interactions (data not shown), and which were hence more likely to have occurred due to linear proximity on the DNA. This resulted in the identification of 6 main interaction classes (Fig. 3a and Extended Data Fig. 7b).

Association of BORIS with lost loops. Only loops that were detected in both the original (sensitive versus resistant) and BORIS depletion (shBORIS versus shCtrl) samples were used for this analysis. First, loops were divided into lost and retained loops upon BORIS depletion, and an odds ratio (two-sided Fisher's exact test) was calculated for the initial presence of BORIS binding on the anchors of these two groups (Fig. 3f). An analogous strategy was followed after first stratifying loops according to the different identified loop classes (Extended Data Fig. 7f, g).

Identification of super-enhancer regions. Super-enhancers were identified using the ROSE algorithm (v.1) (https://bitbucket.org/young_computation/rose). In short, H3K27ac enriched regions were identified with MACS2 and termed typical enhancers. These regions were stitched together if they were within 12.5 kb of each other. Stitched regions were ranked by H3K27ac signal therein and the inclination point at which the two classes of enhancers separated was determined by ROSE. Stitched enhancers above this threshold were considered super-enhancers and the others, typical enhancers. To compare different samples, we used the same maximum threshold between the conditions considered (Extended Data Fig. 8a). Identification of cell-type-specific super-enhancers. Cell-type-specific and active super-enhancers were identified by merging both sensitive- and resistant-cell super-enhancers and determining cell-type specificity based on the differential normalized read density of both H3K27ac and BRD4. In brief, ratios [log2 (resistant/sensitive)] were calculated for H3K27ac and BRD4. A combined threshold of 2.5 was required to identify a cell-type-specific super-enhancer with at least a minimum 0.75 change for each individual mark. Super-enhancers that did not meet these criteria were classed as shared (neutral) between cell types (Extended Data Fig. 8b).

Correlation analysis of looping with gene expression and enhancer landscape. Regulatory interactions were associated to target genes and super-enhancers based on proximity to the TSS and minimal overlap (1 bp) with its anchors, respectively (Fig. 4a and Extended Data Fig. 8f).

Chromatin-based gene classification. Genes were classified as having an 'open', 'neutral' or 'closed' chromatin state based on unsupervised clustering of a metagene representation of ChIP-seq occupancy of H3K27ac and H3K27me3. Each gene (from TSS to TES, and 2 kb up- and downstream of this region) was divided into 20 equally sized bins; the extended regions were binned in regions of 50 bp. Normalized bedgraph files were used to calculate read density (RPM per bp) and *k*-means clustering was applied to group each extended gene region in one of three clusters (Extended Data Fig. 8d, e). An aggregated summary profile was created for each group of genes. The open and closed clusters were classified based on pre-dominantly H3K27ac and H3K27me3 accumulation, respectively, and the 'neutral' cluster displayed on average equal levels of both.

Integrated genomic data analysis. An ensemble analysis was performed to identify the set of genes that showed characteristics of reactivation in resistant cells. For each gene, five features were examined: (i) creation of a unique regulatory interaction; (ii) deposition of BORIS on its promoter or looped enhancer; (iii) association with a resistant cell-specific super-enhancer through overlap with either its promoter or looped anchor; (iv) increased mRNA expression; and (v) transition from a closed or neutral state to an open chromatin state. A unique set of 89 genes (Supplementary Table) that exhibited four out of five features were identified as the top reactivated genes in resistant cells. Within these 89 genes, 13 were identified as transcription factors by the TcoF database (http://www.cbrc. kaust.edu.sa/tcof/) (Fig. 4b).

Allen Brain atlas gene signature. Expression data and metadata for human brain development were downloaded from the Allen Brain atlas (http://www.brainspan. org). Row-normalized *z*-scores of $[log_2(RPKM + 1)]$ values were used to create a heat map. Values greater than 3.5 were set to 3.5 to reduce the effect of extreme outliers on the visualization. Samples were ordered according to developmental time points (Extended Data Fig. 8g).

BORIS and BRD4 correlation at promoter regions. BORIS and BRD4 colocalization and correlation were assessed for the promoter regions of the 89 top-ranked genes. The TSS was extended in both directions by 2 kb and binned in 100-bp regions. Normalized read densities for BORIS and BRD4 were calculated and a Spearman's rank correlation coefficient calculated for sensitive and resistant cells. An aggregated density plot of all 89 genes was created to visualize the increased deposition and correlation of BRD4 and BORIS in resistant cells (Extended Data Fig. 9a).

Gene expression and DNA-binding analysis. To examine the association between gene expression and overlapping targets of MYCN and BORIS in sensitive and resistant cells, respectively, the percentage of gene promoters (± 2 kb TSS) that overlapped with ChIP-seq peaks in 10 equally sized bins based on the expression distribution was calculated (Extended Data Fig. 6f). To visualize and correlate gene expression with DNA binding of MYCN or BORIS, genes were ranked based on expression and plotted against the total rescaled (0–100) binding intensities calculated for each gene promoter (± 2 kb TSS). For each ChIP-seq mark a loess regression curve was computed using a span of 0.1 (Extended Data Fig. 6g).

Transcription factor motif enrichment analysis. Statistically overrepresented motifs were identified with HOMER⁴⁹ (v.2) using the command findMotifs.pl providing both target and background fasta sequences for regions of interest. For promoter regions we selected the top 1,000 up- and downregulated genes in resistant cells and extended the TSS of each gene by 2 kb in both directions. The genomic coordinates were used to extract fasta sequences with the Biostrings package (v.2.50.1) in R and used as target or background to identify motifs associated with promoter regions of genes within each cell type. A similar strategy was followed to identify overrepresented motifs associated with cell-type-specific super-enhancers. Target and background fasta sequences were extracted from the summits of BRD4 peaks located on cell-type-specific super-enhancers and extended by 500 bp in both directions. For a selection of enriched sequences, the associated transcription factor motif and significance level (P) was visualized using a heat map (Fig. 4d). scRNA-seq analysis. The Cell Ranger Single Cell Software Suite, v.1.3 was used to perform sample de-multiplexing, barcode and unique molecular identifier (UMI) processing, and single-cell 3' gene counting. A detailed description of the pipeline and specific instructions to run it can be found at: https://support.10xgenomics. com/single-cell-gene-expression/software/pipelines/latest/what-is-cell-ranger. A high-quality gene expression matrix was created in sequential preprocessing steps. First UMI-based counts were converted to relative expression concentrations by rescaling each cell to a library size of 10,000. Genes were considered detected if rescaled count $> \log_2(0.1 + 1)$ and retained for further analysis if present in at least 0.5% of the cells from the sample with the lowest cell count. Cells were removed if fewer than 1,000 genes were detected. To remove low-quality cells, we calculated five technical indicators (ratio of detected genes/UMI, percentage of mitochondrial genes, percentage of ribosomal genes, average GC content of library and library complexity measured by Shannon Entropy) and performed PCA on indicators with a coefficient of variation > 5%. Next, density-based clustering was performed on the first and second principal component using an epsilon determined by a k-nearest neighbour plot. All cells that were located outside the main cluster were considered low quality and removed from further analysis. Next, we used the R package 'scater' (v.1.10.0) to confirm that there were no technical or experimental confounding effects and the R package 'Seurat' (v.2.3.4) to analyse and visualize the data. In brief, UMI counts were log-normalized with a scale factor of 10,000 and subsequently centre-scaled. To visualize cells in a reduced dimensionality, PCA was performed on the most variable genes, which were identified as genes with higher-than-expected variability in consecutive ranked expression bins. Higher complexity clustering was performed with t-SNE using the first 10 principal components, which were deemed most informative based on heat map and elbow plot observation. To identify homogeneous subpopulations, we performed iterative clustering using the network-based clustering algorithm (shared nearest neighbour) with different resolutions as input until each sample was at least separated in two groups. A simple pseudotime analysis was performed by calculating an average expression profile for each identified subpopulation and ordering them according to the summarized expression of transcription factors that displayed variable expression between sensitive and intermediate or intermediate and resistant cells. Variable expression was defined as showing at least a 33% change in the rank of expression between two samples with a minimal normalized expression level > 0.2. For each sample comparison, at least the top 10 most variable transcription factors were included. In total this resulted in 32 transcription factors. Gene expression values were then linearly rescaled between 0 and 10 to jointly visualize relative expression changes during this pseudotime. To examine co-detection or mutual exclusivity between genes of interest, a two-sided Fisher's exact test was performed for all cells in a given sample. A score combining both the odds ratio and the $-\log_{10}(P \text{ value})$ was calculated to visualize both the strength and direction between genes in pairwise co-expression tests.

Statistical analysis. Analysis for each plot is listed in the figure legend and/or in the corresponding Methods. In brief, all grouped data are presented as mean \pm s.d. unless stated otherwise. All box and whisker plots of expression data are presented as: centre lines, medians; box limits, twenty-fifth and seventy-fifth percentiles; whiskers, minima and maxima (1.5× the interquartile range). Statistical significance for pairwise comparisons was determined using the two-sided Wilcoxon rank-sum test or two-sided unpaired *t*-test, unless stated otherwise. Survival analysis was performed using the Kaplan–Meier method and differences between groups calculated by the two-sided log-rank test and the Bonferroni correction method. Tumour volume comparisons for the xenograft studies were analysed by Mann–Whitney *U* test. *P < 0.05; *P < 0.01. Statistical comparisons of distributions of

fold changes for the expression microarrays were done using the Mann–Whitney U test. All quantitative analyses are expressed as the mean \pm s.d. of three biological replicates, unless stated otherwise. Microarray and ChIP–seq data are based on at least two independent experiments. For all experiments, no statistical methods were used to predetermine sample size. Unless stated otherwise, experiments were not randomized and investigators were not blinded to allocation during experiments and outcome assessment.

Track visualizations. Peaks, (super-) enhancers and HiChIP interactions were visualized with a custom build tool (github.com/RubD/GeTrackViz2) or with the circlize package (v.0.4.5) in R.

Retrospective analysis of gene expression in human samples. Gene expression levels or correlations across primary tumours, healthy tissues or experimental data and patient survival were determined through analysis of the TCGA and TARGET (https://cancergenome.nih.gov/), GTEx (https://www.gtexportal.org/home/), R2 (https://hgserver1.amc.nl/cgi-bin/r2/main.cgi), Allen Brain atlas (http://www.brain-map.org/) and selected datasets representing distinct tumour types with poor prognosis feature annotations (GSE49710 (Neuroblastoma)⁵⁰, GSE17679 (Mixed Ewing Sarcoma)⁵¹, GSE63074 (Non-small cell lung carcinoma)⁵², GSE15709 (ovarian cancer)⁵³, GSE16179 (breast cancer)⁵⁴ and GSE7181 (Glioblastoma)⁵⁵). **Reporting summary.** Further information on research design is available in the Nature Research Reporting Summary linked to this paper.

Data availability

The microarray, ChIP–seq, HiChIP and scRNA-seq datasets generated and analysed during the current study are available in the Gene Expression Omnibus (GEO) repository under accession number GSE103084. The authors declare that all other data supporting the findings of this study are available within the paper and its Supplementary Information files.

Code availability

Custom code is available upon reasonable request.

- Batra, S., Reynolds, C. P. & Maurer, B. J. Fenretinide cytotoxicity for Ewing's sarcoma and primitive neuroectodermal tumor cell lines is decreased by hypoxia and synergistically enhanced by ceramide modulators. *Cancer Res.* 64, 5415–5424 (2004).
- Whang-Peng, J. et al. Cytogenetic characterization of selected small round cell tumors of childhood. *Cancer Genet. Cytogenet.* 21, 185–208 (1986).
- Filippakopoulos, P. et al. Selective inhibition of BET bromodomains. Nature 468, 1067–1073 (2010).
- 34. Marsilje, T. H. et al. Synthesis, structure-activity relationships, and in vivo efficacy of the novel potent and selective anaplastic lymphoma kinase (ALK) inhibitor 5-chloro-N2-(2-isopropoxy-5-methyl-4-(piperidin-4-yl)phenyl)-N4-(2-(isopropylsulfonyl)phenyl)pyrimidine-2,4-diamine (LDK378) currently in phase 1 and phase 2 clinical trials. J. Med. Chem. 56, 5675–5690 (2013).
- Johnson, T. W. et al. Discovery of (10*R*)-7-amino-12-fluoro-2,10,16-trimethyl-15-oxo-10,15,16,17-tetrahydro-2*H*-8,4-(metheno)pyrazolo[4,3-*h*] [2,5,11]-benzoxadiazacyclotetradecine-3-carbonitrile (PF-06463922), a macrocyclic inhibitor of anaplastic lymphoma kinase (ALK) and c-ros oncogene 1 (ROS1) with preclinical brain exposure and broad-spectrum potency against ALK-resistant mutations. *J. Med. Chem.* **57**, 4720–4744 (2014).
- Gosmini, R. et al. The discovery of I-BET726 (GSK1324726A), a potent tetrahydroquinoline ApoA1 up-regulator and selective BET bromodomain inhibitor. J. Med. Chem. 57, 8111–8131 (2014).
- Lovén, J. et al. Revisiting global gene expression analysis. *Cell* 151, 476–482 (2012).
- Debruyne, D. N. et al. ALK inhibitor resistance in ALK(F1174L)-driven neuroblastoma is associated with AXL activation and induction of EMT. Oncogene 35, 3681–3691 (2016).
- Sancak, Y. et al. The Rag GTPases bind Raptor and mediate amino acid signaling to mTORC1. Science 320, 1496–1501 (2008).
- Meerbrey, K. L. et al. The pINDUCER lentiviral toolkit for inducible RNA interference in vitro and in vivo. *Proc. Natl Acad. Sci. USA* **108**, 3665–3670 (2011).
- Tomayko, M. M. & Reynolds, C. P. Determination of subcutaneous tumor size in athymic (nude) mice. Cancer Chemother. Pharmacol. 24, 148–154 (1989).
- Zheng, G. X. et al. Massively parallel digital transcriptional profiling of single cells. *Nat. Commun.* 8, 14049 (2017).
- Gautier, L., Cope, L., Bolstad, B. M. & Irizarry, R. A. affy—analysis of Affymetrix GeneChip data at the probe level. *Bioinformatics* 20, 307–315 (2004).
- Smyth, G. K., Yang, Y. H. & Speed, T. Statistical issues in cDNA microarray data analysis. *Methods Mol. Biol.* 224, 111–136 (2003).

- Bolger, A. M., Lohse, M. & Usadel, B. Trimmomatic: a flexible trimmer for Illumina sequence data. *Bioinformatics* **30**, 2114–2120 (2014).
- Servant, N. et al. HiC-Pro: an optimized and flexible pipeline for Hi-C data processing. *Genome Biol.* 16, 259 (2015).
- Lareau, C. A. & Aryee, M. J. diffloop: a computational framework for identifying and analyzing differential DNA loops from sequencing data. *Bioinformatics* 34, 672–674 (2018).
- Ji, X. et al. 3D chromosome regulatory landscape of human pluripotent cells. Cell Stem Cell 18, 262–275 (2016).
- Heinz, S. et al. Simple combinations of lineage-determining transcription factors prime *cis*-regulatory elements required for macrophage and B cell identities. *Mol. Cell* 38, 576–589 (2010).
- SEQC/MAQC-III Consortium. A comprehensive assessment of RNA-seq accuracy, reproducibility and information content by the Sequencing Quality Control Consortium. Nat. Biotechnol. 32, 903–914 (2014).
- Savola, S. et al. High expression of complement component 5 (C5) at tumor site associates with superior survival in Ewing's sarcoma family of tumour patients. *ISRN Oncol.* 2011, 168712 (2011).
- Huang, S. et al. Analytical performance of a 15-gene prognostic assay for early-stage non-small-cell lung carcinoma using RNA-stabilized tissue. JMD 17, 438–445 (2015).
- Li, M. et al. Integrated analysis of DNA methylation and gene expression reveals specific signaling pathways associated with platinum resistance in ovarian cancer. BMC Med. Genomics 2, 34 (2009).
- Liu, L. et al. Novel mechanism of lapatinib resistance in HER2-positive breast tumor cells: activation of AXL. Cancer Res. 69, 6871–6878 (2009).
- Beier, D. et al. CD133⁺ and CD133⁻ glioblastoma-derived cancer stem cells show differential growth characteristics and molecular profiles. *Cancer Res.* 67, 4010–4015 (2007).

Acknowledgements We thank the George, Young and Gray laboratories and J. R. Gilbert for discussions, and C. Li for assistance with the HiChIP experiments. We thank J. Qi for providing JQ1, and D. Sabatini and S. Elledge for sharing plasmids, pLKO.1 GFP shRNA and plnducer20, respectively. We thank A. Ward and C. Clinton at DFCI Pediatric Oncology and J. Chan and the Clark Smith Tumor Bank, Charbonneau Cancer Institute, Calgary, Canada, for the human tumour samples. We thank Applied Pathology Systems, the DFCI Molecular Biology Core, the Whitehead Genome Technology Core, D. Adeegbe and the NYULH Genome Technology Center for technical support. The results shown here are in part based on data generated by the TCGA Research Network: http://cancergenome.nih.gov/, the R2: Genomics Analysis and Visualization Platform: http://r2.amc.nl/ and the Allan Brain Map Data Portal: http://www.brain-map.org/. This work was supported by NIH grants R01CA197336 (R.E.G. and R.A.Y.), R01CA148688 (R.E.G. and N.S.G.) and a Hyundai Hope on Wheels Scholar Grant (R.E.G.). D.N.D. is a recipient of a Young Investigator Grant from the Alex's Lemonade Stand Foundation/Northwestern Mutual Foundation. D.S.D. is supported by an American Cancer Society fellowship PF-16-146-01-DMC. The NYULH Genome Technology Center is partially supported by the Cancer Center Support Grant P30CA016087 at the Laura and Isaac Perlmutter Cancer Center.


Author contributions D.N.D. and R.E.G. conceived the project and designed the experiments. D.N.D. planned and performed the molecular, cellular and genomic studies. R.D. performed computational analyses with input from D.S.D. and E.M. S.S. contributed to the ChIP-seq and HiChIP experiments. D.S. and S.H.O. contributed to the HiChIP experiment. D.N.D., Y.G. and T.C. performed the mouse experiments. B.S. and M.M. provided technical assistance. H.H. performed the co-immunoprecipitation experiments. L.M. performed the FISH analysis. N.S.G. provided TAE684 and E9. G.-C.Y. supervised the bioinformatics analyses. K.-K.W. enabled the mouse and scRNA-seq studies. D.N.D., R.A., and R.E.G. interpreted the data. D.N.D., R.D. and R.E.G. wrote the manuscript with input from R.A.Y. R.E.G. supervised the research. All authors edited the manuscript.

Competing interests N.S.G. is a founder, SAB member and equity holder of Gatekeeper, Syros Pharmaceuticals, Petra, C4, B2S and Soltego. The Gray laboratory receives or has received research funding from Novartis, Takeda, Astellas, Taiho, Janssen, Kinggen, Voronoi, Her2llc, Deerfield and Sanofi. S.H.O. is a SAB member of Syros. R.A.Y. is a founder and shareholder of Syros, Camp4 Therapeutics, Omega Therapeutics and Dewpoint Therapeutics. R.E.G. is a SAB member of Global Gene Corp.

Additional information

Supplementary information is available for this paper at https://doi. org/10.1038/s41586-019-1472-0.


Correspondence and requests for materials should be addressed to R.E.G. Reprints and permissions information is available at http://www.nature.com/ reprints.



Extended Data Fig. 1 | BORIS is expressed in several cancers and associated with high-risk features. a, b, Relative mRNA expression $[\log_2(FPKM + 1)]$ of CTCF and BORIS in normal tissues (a) and in various cancer types based on TCGA datasets (b). FPKM, fragments per kilobase of transcript per million mapped reads. Keys to cancer types: ACC, adrenocortical carcinoma; AML, acute myeloid leukaemia; BLCA, bladder urothelial carcinoma; BRCA, breast invasive carcinoma; CESC, cervical squamous cell carcinoma and endocervical adenocarcinoma; CHOL, cholangiocarcinoma; COAD, colon adenocarcinoma; DLBC, diffuse large B-cell lymphoma; ESCA, oesophageal carcinoma; GBM, glioblastoma multiforme; HNSC, head and neck squamous cell carcinoma; LGG, low-grade glioma; KICH, kidney chromophobe; KIRC, renal clear cell carcinoma; KIRP, kidney renal papillary cell carcinoma; LAML, acute myeloid leukaemia; LIHC, hepatocellular carcinoma; LUAD, lung adenocarcinoma; LUSC, lung squamous cell carcinoma; MESO, mesothelioma; NB, neuroblastoma; OV, serous ovarian cystadenocarcinoma; PAAD, pancreatic adenocarcinoma; PCPG,


pheochromocytoma and paraganglioma; PRAD, prostate adenocarcinoma; READ, rectum adenocarcinoma; RT, rhabdoid tumour; SARC, sarcoma; SKCM, skin cutaneous melanoma; STAD, stomach adenocarcinoma; TGCT, testicular germ cell tumour; THCA, thyroid carcinoma; THYM, thymoma; UCEC, uterine corpus endometrial carcinoma; UCS, uterine carcinosarcoma; UVM, uveal melanoma; WT, Wilms tumour. c, Box plots showing the correlation of BORIS expression with risk status, tumour stage (primary versus metastasis/recurrence), presence of cancer stem cells (CD133 positivity) and response to targeted (lapatinib) or cytotoxic (cisplatin) therapy in the tumour types depicted. NSCLC, non-small cell lung cancer. Datasets (Mixed Ewing Sarcoma-Savola-117 and NSCLC-Plamadeala-410) were extracted from the R2: Genomics Analysis and Visualization Platform (http://r2.amc.nl). GSE7181 (glioblastoma); GSE16179 (breast cancer); GSE15372 (ovarian cancer). P values determined by two-sided Wilcoxon rank-sum test. For all panels, sample sizes (*n*) are depicted in parenthesis and box plots are as defined in Fig. 4.

RESEARCH LETTER



Extended Data Fig. 2 | ALK inhibitor-resistant cells exhibit stable resistance in vivo and no longer rely on ALK signalling. a, Left, tumour volumes of sensitive and resistant cell xenografts in untreated NU/NU (Crl:NU-Foxn1^{nu}) mice established by subcutaneous injection of 2×10^6 cells into both flanks. Animals were euthanized when tumours reached 1,500–2,000 mm³. Data are mean \pm s.e.m., n = 4 per arm. Right, immunoblot analysis of total and phosphorylated ALK in TAE-resistant xenograft tumours (1 and 2) and sensitive and resistant cells in culture. b, Dose-response curves for TAE684 in sensitive and resistant cell lines established from the same tumour xenografts as in a (IC₅₀ values: sensitive, 7.9 nM; resistant, 878.6 nM). Data are mean \pm s.d., n = 3 biological replicates. c, Tumour volumes (left) and Kaplan-Meier survival curves (right) of resistant cell xenografts in NU/NU (CrTac:NCr-Foxn1^{nu}) mice treated with TAE684 (10 mg kg⁻¹ by oral gavage once daily) or vehicle control for up to 56 days. Data are mean \pm s.e.m., n = 8 per arm. P values determined by Mann–Whitney U test for tumour volumes (P = 0.8404) and by log-rank test for Kaplan–Meier survival analysis (P = 0.8076), both two-sided. d, Dose-response curves for TAE684-sensitive and -resistant

cells treated with ceritinib (IC50 values: sensitive, 33.8 nM; resistant, 446.5 nM) or lorlatinib (IC50 values: sensitive, 47.5 nM; resistant, 2,318 nM). Data are mean \pm s.d., n = 3 biological replicates. **e**, Immunoblot analysis of the indicated proteins in sensitive and resistant cells treated with DMSO or 1 µM TAE684 for 6 or 24 h. f, Electropherograms of ALK kinase domain sequencing in sensitive and resistant cells. Arrows show the F1174L mutation characteristic of Kelly cells. HEK293T cells were used as a control for sequencing wild-type ALK. g, Phosphoproteomic analysis of a panel of receptor tyrosine kinases (RTKs) in sensitive and resistant cells. Each RTK is shown in duplicate and the pairs in the corners of each array are positive controls. Numbered RTKs with corresponding names listed on the right represent the highest-phosphorylated proteins. ALK is depicted in red. h, Quantitative reverse transcription PCR (qRT-PCR) and immunoblot analysis of ABCB1 and ABCG2 multidrug transporter expression in sensitive and resistant cells. The qRT-PCR data are means of n = 2 biological replicates. In **a** (immunoblot), **d**, **f** and **g**, data are representative of two independent experiments (see Supplementary Note 1 for details; for gel source data, see Supplementary Fig. 1).




Extended Data Fig. 3 | See next page for caption.




Extended Data Fig. 3 | Development of resistance is associated with loss of MYCN followed by gradual induction of proneural transcription factors. a, TAE684 dose–response curves of Kelly neuroblastoma cells during resistance establishment (IC₅₀ values: sensitive, 39.4 nM; intermediate, 618 nM; resistant, 1,739 nM). Data are mean \pm s.d., n = 3 biological replicates. Schematic representation of development of resistance is shown above. **b**, *t*-SNE plot of scRNA-seq data showing the segregation of sensitive (n = 5,432), intermediate (n = 6,376) and resistant (n = 6,379) cells. **c**, *t*-SNE plot depicting unsupervised clusters for the individual subpopulations that underlie the pseudotime analysis. **d**, Heat map of rescaled gene expression values of the most variable ranked transcription factors in the three cell states. **e**, qRT–PCR and immunoblot analysis of MYCN expression in TAE684-resistant xenograft tumours (1 and 2) and sensitive and resistant cells in culture. The qRT–PCR data are mean \pm s.d., n = 4 biological replicates for sensitive and resistant cells

(*** $P = 1.396 \times 10^{-11}$; unpaired two-sided *t*-test) and n = 3 technical replicates for each tumour. **f**, Fluorescence in situ hybridization of MYCN in sensitive and resistant cells (representative of 20 nuclei per condition). **g**, ChIP-seq track of H3K27me3 binding at the *MYCN* locus in sensitive and resistant cells. Signal intensity is given in the top right corner. **h**, Line plot showing the association between genes ordered by expression (*x* axis) and changes in absolute gene expression levels (*y* axis) in sensitive or resistant cells. Bar plot, total transcriptional yield in sensitive or resistant cells. **i**, Immunoblot analysis of the indicated proteins in sensitive and resistant cells expressing control (shCtrl) or MYCN (shMYCN-1 and -2) shRNAs. **j**, Violin plots representing the expression distribution of selected genes in the same cells as in **a** (centre line, median). **k**, Bar plot showing the fractions of cells with detectable mRNA levels of the same genes as in **d**. In **e** (immunoblot) and **f**-**i**, data are representative of two independent experiments (for gel source data, see Supplementary Fig. 1).



Extended Data Fig. 4 | Overexpression of *BORIS* is seen in resistance models of neuroblastoma and correlates with high-risk disease and a poor outcome. a, qRT–PCR and immunoblot analysis of BORIS expression in TAE684-resistant Kelly cell xenograft tumours (1 and 2) and sensitive and resistant cells in culture. The qRT–PCR data are mean \pm s.d., n = 4 biological replicates for sensitive and resistant cells (**P = 0.0014; unpaired two-sided *t*-test) and n = 3 technical replicates for each tumour. b, Bisulfite sequencing of the *BORIS* promoter in sensitive and resistant cells. Black circles represent methylated cytosines. The B and C TSSs are indicated by arrows. c, Dose–response curves to TAE684 (left) and immunoblot analysis of BORIS expression (right) in TAE684-sensitive and -resistant SK-N-SH neuroblastoma cells (IC₅₀ values: sensitive, 47.9 nM; resistant, 1,739 nM). d, Dose–response curves to the CDK12 inhibitor, E9 (left) and immunoblot analysis of BORIS expression (right) in sensitive

and resistant SK-N-BE(2) neuroblastoma cells (IC₅₀ values: sensitive, 9.5 nM; resistant, 638 nM). Data are mean \pm s.d., n = 3 biological replicates for **c** (left) and **d** (left). **e**, Immunohistochemical staining of BORIS expression in primary neuroblastoma tumour samples (representative of four independent experiments). Scale bar, 20 µm. **f**, Box plots showing correlation of BORIS expression with the indicated parameters in a human neuroblastoma dataset (n = 498; Tumour Neuroblastoma-SEQC-498; R2: Genomics Analysis and Visualization Platform (http://r2.amc.nl)). Box plots are as defined in Fig. 4. *P* values were determined by two-sided Wilcoxon rank-sum test. **g**, Kaplan-Meier analysis of overall survival based on BORIS expression in the same dataset as in **f** (n = 498; two-sided log-rank test with Bonferroni correction). In **a**, **c**, **d** (immunoblots) and **b**, data are representative of two independent experiments. Sample sizes (n) are depicted in parenthesis for **f** and **g** (for gel source data, see Supplementary Fig. 1).




Extended Data Fig. 5 | Resistant cells are dependent on BORIS for survival. a, Dose-response curves to TAE684 in resistant cells expressing control (shCtrl) or BORIS (shBORIS) shRNAs (IC₅₀ values: shCtrl, 537.7 nM; shBORIS, 141.2 nM). Data are mean \pm s.d., n = 3 biological replicates. b, Heat map of gene expression values in the same cells as in a (n = 2 biological replicates). Rows are *z*-scores calculated for each gene in both conditions. c, Immunoblot analysis of the indicated proteins in the same cells as in a. d, e, Immunoblot analysis of the indicated proteins (Cl., cleaved; CC3, cleaved caspase 3) (d), and quantification of trypan blue staining (e) in sensitive and resistant cells expressing control (shCtrl) or BORIS (shBORIS-3 and -4) shRNAs. Data are mean \pm s.d., n = 3

RESEARCH

LETTER


biological replicates (*P < 0.05; **P < 0.01; ***P < 0.001; unpaired two-sided *t*-tests). **f**-**h**, Phase-contrast microscopy images (scale bars, 150 µm) (**f**), growth curves (**g**) and flow cytometry analyses (**h**) of propidium iodide (PI) staining in sensitive, intermediate and resistant cells. Data are mean \pm s.d., n = 3 biological replicates (***P < 0.0001 for all comparisons; two-way ANOVA). **i**, qRT–PCR analysis of the expression of the indicated proneural transcription factors in the same sensitive (DMSO) versus MYCN^{KD} and BORIS^{Ind} (DOX + TAE) cells as in Fig. 1g. Data are mean \pm SD, n = 3 biological replicates (*P < 0.05; **P < 0.01; unpaired two-sided *t*-tests). In **c**, **d**, **f** and **h**, data are representative of two independent experiments (for gel source data, see Supplementary Fig. 1).








Extended Data Fig. 6 | BORIS colocalizes with CTCF and open chromatin. a, Bar graphs illustrating the overlap of shared and specific BORIS and CTCF-binding sites in sensitive and resistant cells. Most resistant cell-specific BORIS peaks (red) colocalize with CTCF peaks that are shared between the two cell types. The markedly lower number of BORIS peaks that are unique to sensitive cells (green) or shared between sensitive and resistant cells (grey) typically do not overlap with CTCF peaks that are shared or specific to any cell type (top). Most CTCF peaks are shared (grey) between sensitive and resistant cells and either do not overlap with BORIS peaks, or overlap only with those restricted to resistant cells (bottom). b, Comparison of CTCF and BORIS peaks identified in sensitive and resistant cells. c, Co-immunoprecipitation of BORIS with CTCF in sensitive and resistant cells (representative of two independent experiments). IgG and sample without antibody (Ab) serve as controls. d, Pie charts depicting the percentages of genomic regions bound by BORIS in sensitive (top) and resistant (bottom) cells. Numbers of BORIS-binding peaks in each cell type are given below each pie chart. The regions shown are promoters (TSS \pm 2 kb), typical enhancers (H3K27ac), active enhancers (H3K27ac + BRD4), repressed chromatin (H3K27me3), exons, introns, and other (peaks not assigned to any of the previous categories). **e**, Meta-analysis of average CTCF and BORIS ChIP-seq signals in RPM per bp at enhancer and TSS regions in sensitive and resistant cells. **f**, Percentage of gene promoters bound by BORIS in sensitive (black) and resistant (red) cells for 10 equal-sized groups ordered based on absolute gene expression levels in resistant cells. Percentage of promoters bound by BORIS in resistant cells that were also originally bound by MYCN in sensitive cells is shown in grey. **g**, Loess regression analysis of ranked gene expression against BORIS and MYCN occupancies at gene promoters in sensitive and resistant cells. Shaded regions represent 95% confidence intervals. All panels except **c** depict data from n = 2biological replicates (for gel source data, see Supplementary Fig. 1).



Extended Data Fig. 7 | See next page for caption.



Extended Data Fig. 7 | Regulatory loops in resistant cells are more vulnerable to BORIS depletion. a, Heat map depicting the Spearman correlation between HiChIP biological replicates of sensitive and resistant cells in genome-wide bins of 5 kb for all merged anchor regions. b, Box plots showing the genomic length distribution (in log₂(bp)) for interaction classes that are specific to resistant cells. c, Table depicting HiChIP loop class statistics in resistant cells, including their association with BORIS binding. d, ChIP-seq tracks of the indicated proteins in sensitive and resistant cells at the *TCP11L2* locus (representative of two independent experiments), with resistant cell-specific regulatory interactions shown below (HiChIP Res: PET numbers, next to each interaction). Signal intensity is given in the top left corner for each track. e, Heat map depicting the Spearman correlation between HiChIP biological replicates of sensitive, resistant, shCtrl and shBORIS cells in genome-wide bins of 5 kb for all merged anchor regions. f, Bar plots showing the number and fraction of resistant cell-specific loops for all interaction classes that were BORIS negative and positive in resistant cells, and that were lost after BORIS depletion. **g**, Bar plots showing the odds ratio (two-sided Fisher's exact test) of losing loops that were previously bound by BORIS for all interaction classes. **h**, Box plots showing the initial intensities (in normalized read counts) of BORIS and SMC1A binding in the shRNA control cells at the anchors of the resistant cell-specific loops that were significantly lost versus those that were retained in shBORIS cells (two-sided Wilcoxon rank-sum test). **i**, Box plot showing the difference in SMC1A loss (shBORIS versus shCtrl) on the same anchors as in **h**. *P* value determined by two-sided Wilcoxon rank-sum test. All box plots are as defined in Fig. 4. **j**, Metaplots depicting BORIS, SMC1A and CTCF binding at the anchors of the resistant cell-specific loops that were lost or retained after BORIS depletion. In **a-c** and **e-g**, data are from n = 3biological replicates. In **h-j**, data are from n = 2 biological replicates.




Extended Data Fig. 8 | See next page for caption.



Extended Data Fig. 8 | Redistribution of the super-enhancer landscape with subsequent expression of a BORIS-dependent proneural network in resistant cells. a, Accumulation of H3K27ac signal at enhancer regions. Typical enhancers (grey) are plotted according to increasing levels of normalized H3K27ac signal (length \times density) in sensitive and resistant cells. The highest cut-off based on the inclination point in both sensitive and resistant cells was used to delineate super-enhancers (red). b, Scatter plot showing differential binding of H3K27ac [log₂(RPM per bp (+ 1) and BRD4 [log₂(RPM per bp + 1)] for all detected super-enhancers in both sensitive and resistant cells. Cell-specific super-enhancers were identified based on the combined increase in H3K27ac and BRD4 binding. For each individual histone mark, a 0.75 log₂-transformed fold change threshold was applied and a minimum summed 2.5 log₂transformed fold change was used as the final cut-off. c, Bar plot depicting the enrichment (two-sided Fisher's exact test) and fractions of resistant cell-specific and shared super-enhancers that were located at resistant cell-specific regulatory loop anchors in resistant cells. d, Density plots showing the aggregated accumulation of H3K27ac and H3K27me3 at gene regions, defined as 2 kb upstream of the TSS and 2 kb downstream of the transcription end site (TES). k-means clustering (k = 3) analysis resulted in the separation of genes associated with 'open', 'neutral' or 'closed' chromatin in both sensitive and resistant cells. e, Sankey diagram of the

distribution of genes in distinct chromatin states and the switches between sensitive and resistant cells. f, Box plots showing the expression level changes upon BORIS depletion for genes that had a resistant cell-specific and BORIS-positive regulatory interaction and were not associated with a super-enhancer (n = 720), associated with a super-enhancer in both cell types (n = 514) or associated with a super-enhancer seen only in the resistant cells (n = 134) (two-sided Wilcoxon rank-sum test). Box plots are as defined in Fig. 4. g, Heat map of the expression levels of the indicated proneural transcription factor genes during brain development (http://www.brain-map.org/). Gene expression levels are represented as *z*-scores for different developmental time points (n = 413; pcw, postconceptional weeks). h, Heat map showing the odds ratios (two-sided Fisher's exact test) for co-detection of the indicated transcription factors based on the scRNA-seq data in resistant cells (n = 6,379). i, Immunoblot analysis of the indicated proteins in sensitive and resistant cells expressing control (shCtrl) or BORIS (shBORIS-3 and -4) shRNAs. j, k, qRT-PCR analysis of the indicated genes (j) and ChIP-qPCR analysis of BORIS binding at the promoter regions of BORIS and NEUROG2 (k) in sensitive and resistant SK-N-BE(2) neuroblastoma cells. Data are mean \pm s.d., n = 3 biological replicates in **j** and **k** (*P < 0.05; **P < 0.01; ***P < 0.001; unpaired two-sided *t*-tests). All other panels except **g** and **h** depict data from n = 2 biological replicates (for gel source data, see Supplementary Fig. 1).




Extended Data Fig. 9 | See next page for caption.



Extended Data Fig. 9 | The proneural transcription factor network in resistant cells is sensitive to BRD4 inhibition. a, Metaplots showing the correlation between BRD4 and BORIS co-occupancies at the promoter regions (\pm 2 kb) of the 89 top-ranked genes in resistant versus sensitive cells based on the features in Fig. 4b (r, Spearman correlation coefficient). b, Immunoblot analysis of BRD4 and cleaved PARP expression in sensitive and resistant cells expressing control (shCtrl) or BRD4 (shBRD4-A and -B) shRNAs. c, Immunoblot analysis of the indicated proteins in sensitive and resistant cells treated with DMSO, TAE684 (1 μ M) or JQ1 (2.5 μ M) for 48 h. d, Dose-response curves for sensitive and resistant cells treated with JQ1 or I-BET726 (JQ1 (IC50 values: sensitive, 4,798 nM; resistant, 645 nM); I-BET726 (IC50 values: sensitive, 6,203 nM; resistant, 347 nM)). Data are mean \pm s.d., n = 3 biological replicates. **e**, Box plots comparing the expression of the transcription factors listed in Fig. 4b (n = 13) with that of all genes (n = 18,038) in sensitive versus resistant cells (left), and between DMSO and JQ1-treated resistant cells (right) (P values determined by two-sided Wilcoxon rank-sum test). f, ChIP-seq tracks of the indicated proteins at the SIX1 or SIX4 locus in sensitive, resistant and JQ1-treated resistant cells (2.5 µM for 48 h). Super-enhancers are depicted as coloured rectangles below the tracks. Signal intensity is shown in the top left corner for each track. g, h, Tumour volumes (g) and survival curves

(h) in sensitive- and resistant-cell xenografts in NU/NU (Crl:NU-Foxn1^{nu}) mice treated with JQ1 (50 mg kg⁻¹ i.p. once daily) and vehicle control for up to 87 days. Data are mean \pm s.e.m., n = 6 per arm. Significance was calculated by Mann–Whitney U test for tumour volumes (sensitive: P = 0.3231; resistant: P = 0.0023) and by log-rank test for Kaplan–Meier survival analysis (sensitive: P = 0.3047; resistant: 0.0348), both two-sided. i, Heat map of gene expression values in sensitive, resistant and JQ1-treated resistant cells. Rows are z-scores calculated for each gene in each condition. j, Number of transcripts in sensitive, JQ1-treated resistant, shBORIS-expressing resistant and resistant cells based on expression array data after spike-in normalization. k, Scatter plot displaying the medianscaled fold-change gene expression values for shBORIS and JQ1-treated resistant cells. The top-ranked transcription factors that show decreased expression levels after both BORIS knockdown and JQ1 treatment are listed in red (bottom left quadrant). The pie chart represents the fraction of all top-ranked transcription factors that are located in the left lower quadrant of the scatter plot. All box plots are as defined in Fig. 4. In b, c and f, data are representative of two independent experiments. In a, e and **i**–**k**, data are from n = 2 biological replicates (see Supplementary Note 2 for further details; for gel source data, see Supplementary Fig. 1).



Extended Data Fig. 10 | Aberrantly expressed BORIS binds to regulatory regions and is associated with new super-enhancers in Ewing

sarcoma cells. a, Immunoblot analysis of BORIS expression in TC-32 (pre-chemotherapy), TC-71 and CHLA-10 (relapsed, post-chemotherapy) Ewing sarcoma cells, compared with BORIS expression in resistant (Kelly) neuroblastoma cells. b, Meta-analysis of average BORIS ChIP-seq signals in RPM per bp at all combined BORIS-binding sites for TC-32 and TC-71 cells. c, Meta-analysis of average BORIS, H3K27ac and SMC1A ChIP-seq signals in RPM per bp at TC-71-specific BORIS-binding sites. d, Pie chart depicting the proportions of genomic regions bound by BORIS in TC-71 cells. The regions shown are promoters (TSS \pm 2 kb), typical and superenhancers (H3K27ac), and other (if peaks were not assigned to any of the previous categories). **e**, Bar plot showing the odds ratios (two-sided Fisher's exact test) of BORIS localization to regulatory genomic regions in TC-71 cells. All panels are representative of two independent experiments (see Supplementary Note 3 for further details; for gel source data, see Supplementary Fig. 1).

natureresearch

Corresponding author(s): Rani E. George

Last updated by author(s): Jun 13, 2019

Reporting Summary

Nature Research wishes to improve the reproducibility of the work that we publish. This form provides structure for consistency and transparency in reporting. For further information on Nature Research policies, see <u>Authors & Referees</u> and the <u>Editorial Policy Checklist</u>.

Statistics

For	all st	atistical analyses, confirm that the following items are present in the figure legend, table legend, main text, or Methods section.
n/a	Cor	firmed
	\square	The exact sample size (n) for each experimental group/condition, given as a discrete number and unit of measurement
	\boxtimes	A statement on whether measurements were taken from distinct samples or whether the same sample was measured repeatedly
		The statistical test(s) used AND whether they are one- or two-sided Only common tests should be described solely by name; describe more complex techniques in the Methods section.
\boxtimes		A description of all covariates tested
	\square	A description of any assumptions or corrections, such as tests of normality and adjustment for multiple comparisons
		A full description of the statistical parameters including central tendency (e.g. means) or other basic estimates (e.g. regression coefficient) AND variation (e.g. standard deviation) or associated estimates of uncertainty (e.g. confidence intervals)
	\boxtimes	For null hypothesis testing, the test statistic (e.g. <i>F</i> , <i>t</i> , <i>r</i>) with confidence intervals, effect sizes, degrees of freedom and <i>P</i> value noted Give <i>P</i> values as exact values whenever suitable.
\boxtimes		For Bayesian analysis, information on the choice of priors and Markov chain Monte Carlo settings
\boxtimes		For hierarchical and complex designs, identification of the appropriate level for tests and full reporting of outcomes
	\square	Estimates of effect sizes (e.g. Cohen's d, Pearson's r), indicating how they were calculated
	,	Our web collection on <u>statistics for biologists</u> contains articles on many of the points above.

Software and code

Policy information about availability of computer code

Data collection	Data collection was done with BD CellQuest Pro software (BD Biosciences) for the flow cytometry experiment. The Cell Ranger Single Cell Software Suite (v1.3) was used to perform sample de-multiplexing, barcode and UMI processing, and singlecell 3' gene counting for the single cell RNA-sequencing experiment.
Data analysis	R (v3.1.3), affy (v1.44.0), limma (v3.22.7), Fastqc (v0.11.5), STAR (v2.5.1b_modified), Samtools (v1.3.1), MarkDuplicates (v2.1.1), deepTools (v2.2.4), MACS2 (v2.1.1), bedGraphToBigWig (v4), R (v3.5.1), Rstudio (v1.1.463), data.table (v1.12.2), trimmomatic (v0.36), HiC-Pro (v2.10.0), hichipper (v0.7.3), diffloop (v1.10.0), ROSE (v1), HOMER (v2) (http://homer.ucsd.edu/homer/motif/), Biostrings (v2.50.1),Cell Ranger (v1.3), Scater (v1.10.0), Seurat (v2.3.4), Circlize (v0.4.5), IGV (v2.3.74), FlowJo (v10.0.5) and GraphPad Prism (v7.02). Custom code is available upon reasonable request.

For manuscripts utilizing custom algorithms or software that are central to the research but not yet described in published literature, software must be made available to editors/reviewers. We strongly encourage code deposition in a community repository (e.g. GitHub). See the Nature Research guidelines for submitting code & software for further information.

Data

Policy information about availability of data

All manuscripts must include a data availability statement. This statement should provide the following information, where applicable:

- Accession codes, unique identifiers, or web links for publicly available datasets
- A list of figures that have associated raw data
- A description of any restrictions on data availability

The microarray, ChIP-seq, HiChIP and scRNA-seq datasets generated and analyzed during the current study are available in the Gene Expression Omnibus repository under accession number GSE103084. The authors declare that all other data supporting the findings of this study are available within the paper and its Supplementary Information files.

Field-specific reporting

Life sciences

Please select the one below that is the best fit for your research. If you are not sure, read the appropriate sections before making your selection.

🗌 Behavioural & social sciences 👘 Ecological, evolutionary & environmental sciences

For a reference copy of the document with all sections, see nature.com/documents/nr-reporting-summary-flat.pdf

Life sciences study design

All studies must disclose on these points even when the disclosure is negative.

Sample size	No statistical methods were used to predetermine sample sizes. Sample sizes were chosen in order to be able to perform statistical analyses, as is standard in the field
Data exclusions	No data were excluded from the analyses.
Replication	To verify the reproducibility of our findings, experiments were performed using at least three biological replicates, unless clearly stated otherwise in the figure legends. All attempts at replication were successful.
Randomization	For in vivo mice experiments, mice were randomized into groups of equal average volumes.
Blinding	For experiments involving human research participants, investigators were blinded to group allocation during data collection. For in vivo mice experiments, investigators were not blinded to group allocation during data collection. For experiments involving cell culture, blinding did not apply.

Reporting for specific materials, systems and methods

We require information from authors about some types of materials, experimental systems and methods used in many studies. Here, indicate whether each material, system or method listed is relevant to your study. If you are not sure if a list item applies to your research, read the appropriate section before selecting a response.

Materials & experimental systems

n/a	Inv	olved in the study
	\boxtimes	Antibodies
	\boxtimes	Eukaryotic cell lines
\boxtimes		Palaeontology
	\boxtimes	Animals and other organisms
	\boxtimes	Human research participants
\boxtimes		Clinical data

Antibodies

Antibodies used The following antibodies were used (all commercially available): From Cell Signaling Technology (Danvers, MA, USA): 1. N-Myc (#9405, lot 2, 1:1,000); 2. N-Myc (#51705, Clone D4B2Y, lot 1); 3. Cleaved PARP Asp214 (#9541, lot 15, 1:4,000); 4. Cleaved caspase 3 Asp175 (#9661, lot 45, 1:500); 5. ALK (#3333, Clone C26G7, lot 7, 1:1,000); 6. Akt (pan) (#4691. Clone C67E7. lot 20. 1:4.000): 7. Phospho-Akt Thr308 (#9275, lot 26, 1:1,000); 8. Phospho-Akt Ser473 (#9271, lot 14, 1:2,000); 9. p44/42 MAPK (ERK1/2) (#4695, Clone 137F5, lot 21, 1:4,000); 10. Phospho- p44/42 MAPK (ERK1/2) (Thr202/Tyr204) (#4377, Clone 197G2, lot 10, 1:4,000); 11. S6 Ribosomal Protein (#2217, Clone 5G10, lot 5, 1:8,000); 12. Phospho-S6 Ribosomal Protein (Ser235/236) (#4857, Clone 91B2, lot 2, 1:8,000); 13. Stat3 (#4904, Clone 79D7, lot 7, 1:2,000); 14. Phospho-Stat3 (Tyr705) (#9131, lot 30, 1:1,000); 15. MDR1/ABCB1 (#12683, Clone D3H1Q, lot 2, 1:5,000); 16. Sox2 (#3579, Clone D6D9, lot 8, 1:2,000); 17. β-Actin (#4967, multiple lots, 1:5,000); 18. CTCF (#3417, Clone D1A7, lot 1); 19. Normal Rabbit IgG (#2729, lot 7)

Methods

n/a	Involved in the study
	ChIP-seq
\square	Flow cytometry

MRI-based neuroimaging

20. Anti-mouse IgG, HRP-linked (#7	076, lot 33, 1:2,000).
From Santa Cruz Biotechnology (San	
21. Anti-rabbit IgG-HRP (sc-2357, m	nultiple lots, 1:5,000).
From Bethyl Laboratories (Montgor	
22. BRD4 (A301-985A100, lot 6, 1:3	10,000);
23. SMC1A (A300-055A, lot 6).	
From Millipore (Billerica, MA, USA):	
24. CTCF (#07-729, lot 2887267, 1:5	
25. SOX9 (#AB5535, lot 2847051, 1	
26. H3K27me3 (#07-449, lot 29728	64).
From Abcam (Cambridge, MA, USA)	
27. ALK (phospho Y1507) (ab73996	
28. BORIS (ab187163, Clone EP1220	
29. H3K27ac (ab2729, lot GR31988)	66-1).
From NOVUS Biologicals (Littleton,	
30. BORIS (NBP2-52405, Clone 2083	11, lot CRT/17/124).
From Active Motif (Carlsbad, CA, US	
31. BORIS (#39851, lot 18916002, 1	1:3,000).
From Sigma-Aldrich (Saint Louis, Mo	
32. SIX1 (HPA001893, lot I114662, 2	1:4,000).
From Abbott (Abbott Park, IL, USA):	
33. Vysis LSI N-MYC (2p24) Spectrui	mGreen/Vysis CEP 2 SpectrumOrange Probe (07J72-001).
,	le antibodies can be found using the links below:
1. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc 5. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc 5. https://www.cellsignal.com/proc 6. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc 5. https://www.cellsignal.com/proc 6. https://www.cellsignal.com/proc 7. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 4. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 4. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 4. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 4. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod mab/4377	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 bducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit-
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 4. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod mab/4377 11. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- bucts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc 5. https://www.cellsignal.com/proc 6. https://www.cellsignal.com/proc 8. https://www.cellsignal.com/proc 9. https://www.cellsignal.com/proc 10. https://www.cellsignal.com/proc mab/4377 11. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/3333 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p4-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 boducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 4. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod mab/4377 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/phospho-akt-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- pducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/4595 oducts/primary-antibodies/phospho-p64-67-rabbit-mab/2217
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- pducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 ducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc 5. https://www.cellsignal.com/proc 6. https://www.cellsignal.com/proc 8. https://www.cellsignal.com/proc 9. https://www.cellsignal.com/proc 10. https://www.cellsignal.com/proc 10. https://www.cellsignal.com/proc 11. https://www.cellsignal.com/proc 12. https://www.cellsignal.com/proc 13. https://www.cellsignal.com/proc 14. https://www.cellsignal.com/proc	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 oducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-erk12-137f5-rabbit-mab/4695 oducts/primary-antibodies/p44-42-mapk-erk12-137f5-rabbit-mab/4695 oducts/primary-antibodies/p44-42-mapk-erk12-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 oducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/mospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- bducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 bducts/primary-antibodies/stat3-79d7-rabbit-mab/4904 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bhospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bhospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bhospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bhospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bhospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bhospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/bactin-antibody/4967
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 18. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/p6-pibospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- bducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 bducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/s02-d6d9-xp-rabbit-mab/12683 bducts/primary-antibodies/s02-d6d9-xp-rabbit-mab/3579 bducts/primary-antibodies/b-actin-antibody/4967 bducts/primary-antibodies/b-actin-antibody/4967
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 oducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 oducts/primary-antibodies/s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/12683 oducts/primary-antibodies/ox2-d6d9-xp-rabbit-mab/3579 oducts/primary-antibodies/b-actin-antibody/4967 oducts/primary-antibodies/b-actin-antibody/4967 oducts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3417 oducts/primary-antibodies/normal-rabbit-igg/2729
1. https://www.cellsignal.com/proc 2. https://www.cellsignal.com/proc 3. https://www.cellsignal.com/proc 4. https://www.cellsignal.com/proc 5. https://www.cellsignal.com/proc 6. https://www.cellsignal.com/proc 7. https://www.cellsignal.com/proc 9. https://www.cellsignal.com/proc 10. https://www.cellsignal.com/proc 12. https://www.cellsignal.com/proc 13. https://www.cellsignal.com/proc 13. https://www.cellsignal.com/proc 15. https://www.cellsignal.com/proc 16. https://www.cellsignal.com/proc 15. https://www.cellsignal.com/proc 16. https://www.cellsignal.com/proc 17. https://www.cellsignal.com/proc 18. https://www.cellsignal.com/proc 19. https://www.cellsignal.com/proc 19. https://www.cellsignal.com/proc 19. https://www.cellsignal.com/proc 10. https://www.cellsignal.com/proc 10	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e-rabbit-mab/3333 ducts/primary-antibodies/akt-pan-c67e-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/4695 oducts/primary-antibodies/s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 bducts/primary-antibodies/stat3-79d7-rabbit-mab/4904 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/12683 bducts/primary-antibodies/ox2-d6d9-xp-rabbit-mab/3579 bducts/primary-antibodies/ox2-d6d9-xp-rabbit-mab/3417 bducts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3417 bducts/primary-antibodies/normal-rabbit-igg/2729 bducts/primary-antibodies/normal-rabbit-igg/2729
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 18. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/3333 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 boducts/primary-antibodies/s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 boducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 boducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 boducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3579 bducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3417 bducts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3417 bducts/primary-antibodies/normal-rabbit-igg/2729 bducts/primary-antibodies/normal-rabbit-igg/2729
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 bducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 bducts/primary-antibodies/stat3-79d7-rabbit-mab/4904 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 bducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/12683 bducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3579 bducts/primary-antibodies/ctcf-d1a7-xp-rabbit-mab/3417 bducts/primary-antibodies/normal-rabbit-igg/2729 bducts/primary-antibodies/normal-rabbit-igg/2729 bducts/secondary-antibodies/normal-rabbit-igg-hrp-linked-antibody/7076 bducts/perimary-antibodies/normal-rabbit-igg-hrp-linked-antibody/7076
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 18. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 oducts/primary-antibodies/s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/so2-d6d9-xp-rabbit-mab/2579 oducts/primary-antibodies/bospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/bospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/bospho-attin-antibody/4967 oducts/primary-antibodies/ctcf-d1a7-xp-rabbit-mab/3417 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/secondary-antibodies/normal-rabbit-igg/2729 oducts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 oducts/psAJ100/BRD4+Antibody ict/A301-985A100/BRD4+Antibody
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 8. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10	Aucts/primary-antibodies/n-myc-antibody/9405 Aucts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Aucts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 Aucts/primary-antibodies/akt-pan-c67e7-rabbit-mab/4691 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9275 Aucts/primary-antibodies/phospho-akt-erk1-2-137f5-rabbit-mab/4695 Aucts/primary-antibodies/phospho-akt-erk1-2-137f5-rabbit-mab/4695 Aucts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 aducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/2683 aducts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3579 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/secondary-antibodies/normal-rabbit-igg/2729 aducts/secondary-antibodies/nati-mouse-igg-hrp-linked-antibody/7076 aducts/secondary-antibodies/nati-mouse-igg-hrp-linked-antibody/7076 aducts/secondary-antibodies/nati-mouse-igg-hrp-linked-antibody/7076 aducts/secondary-antibodies/nati-mouse-igg-hrp-linked-antibody/7076 aducts/secondary-antibodies/nati-mouse-igg-hrp-linked-antibody/7076 aducts/secondary-antibodies/nati-mouse-igg-hrp-linked-antibody/7076 aducts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/a00-055A/SMC1+Antibody (US/en/product/Anti-CTCF-Antibody,MM_NF-07-729
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 18. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 1	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/phospho-akt-thr308-antibody/9271 ducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 bducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 bducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 bducts/primary-antibodies/sta3-79d7-rabbit-mab/4904 bducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/4904 bducts/primary-antibodies/cfcf-d1a7-xp-rabbit-mab/3579 bducts/primary-antibodies/cfcf-d1a7-xp-rabbit-mab/3579 bducts/primary-antibodies/normal-rabbit-igg/2729 bducts/primary-antibodies/normal-rabbit-igg/2729 bducts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 bducts/primary-antibodies/normal-rabbit-gg-2729 bducts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 bducts/primary-antibodies/normal-rabbit-gg/2729 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/primary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/secondary-antibodies/notf-cf-d1a7-xp-rabbit-mab/3579 bducts/s
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 18. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 21. https://www.cellsignal.com/prod 21. https://www.cellsignal.com/prod 21. https://www.cellsignal.com/prod 21. https://www.cellsignal.com/prod 22. https://www.cellsignal.com/prod 23. https://www.cellsignal.com/prod 24. https://www.cellsignal.com/prod 25. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 27. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 27. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 27. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 26. https://www.cellsignal.com/prod 27. https://www.cellsignal.com/prod 2	ducts/primary-antibodies/n-myc-antibody/9405 ducts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 ducts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 ducts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 ducts/primary-antibodies/alk-c26g7-rabbit-mab/3333 ducts/primary-antibodies/alk-c26g7-rabbit-mab/4691 ducts/primary-antibodies/phospho-akt-thr308-antibody/9275 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-akt-ser473-antibody/9271 ducts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- boducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/s6-ribosomal-protein-5g10-rabbit-mab/2217 oducts/primary-antibodies/s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/s6-ribosomal-protein-sg10-rabbit-mab/2217 oducts/primary-antibodies/sospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 oducts/primary-antibodies/sospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/12683 oducts/primary-antibodies/b-actin-antibody/4967 oducts/primary-antibodies/ctcf-d1a7-xp-rabbit-mab/3579 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/secondary-antibodies/normal-rabbit-igg/2729 oducts/secondary-antibodies/normal-rabbit-igg/2729 oducts/secondary-antibodies/normal-rabbit-igg/2729 vducts/secondary-antibodies/normal-rabbit-igg/2729 vducts/secondary-antibodies/normal-rabbit-igg/2729 vducts/secondary-antibodies/normal-rabbit-igg/2729 vducts/secondary-antibodies/normal-rabbit-igg/2729 vducts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 oduct/nouse-anti-rabbit-igg-hrp vct/A301-985A100/BRD4+Antibody v(Js/en/product/Anti-CTCF-Antibody,MM_NF-07-729 v(Js/en/product/Anti-Sox9-Antibody,MM_NF-07-729 v(Js/en/product/Anti-Sox9-Antibody,MM_NF-07-7449
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 20. https://www.cellsignal.com/prod 21. https://www.cellsignal.com/prod 23. https://www.cellsignal.com/prod 24. http://www.endmillipore.com/ 25. http://www.abcam.com/alk-pf	Aucts/primary-antibodies/n-myc-antibody/9405 Aucts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Aucts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 Aucts/primary-antibodies/alk-c26g7-rabbit-mab/4691 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9275 Aucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Aucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Aucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Aucts/primary-antibodies/phospho-p44-42-mapk-erk1-2-137f5-rabbit-mab/4695 aducts/primary-antibodies/phospho-set-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-set-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-set-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-set-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/stat3-79d7-rabbit-mab/4904 aducts/primary-antibodies/stat3-79d7-rabbit-mab/4904 aducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/12683 aducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3579 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-trimethyl-Histone-H3-Lys27-Antibody,MM_NF-07-449 hospho-y1507-antibody-ab73996.html
1. https://www.cellsignal.com/prod 2. https://www.cellsignal.com/prod 3. https://www.cellsignal.com/prod 5. https://www.cellsignal.com/prod 6. https://www.cellsignal.com/prod 7. https://www.cellsignal.com/prod 9. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 10. https://www.cellsignal.com/prod 11. https://www.cellsignal.com/prod 12. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 13. https://www.cellsignal.com/prod 14. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 15. https://www.cellsignal.com/prod 16. https://www.cellsignal.com/prod 17. https://www.cellsignal.com/prod 19. https://www.cellsignal.com/prod 20. https://www.cellsignal.com/prod 21. https://www.cellsignal.com/prod 22. https://www.cellsignal.com/prod 23. https://www.seta.com/soth/prod 24. http://www.endmillipore.com/ 25. http://www.seta.com/soth/soth/ 25. https://www.seta.com/soth/soth/ 26. https://www.seta.com/soth/soth/soth/ 27. https://www.seta.com/soth/soth/soth/ 28. https://www.seta.com/soth/soth/soth/soth/soth/ 28. https://www.seta.com/soth/soth/soth/soth/soth/soth/soth/soth	Aucts/primary-antibodies/n-myc-attibody/9405 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Aucts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 Aucts/primary-antibodies/alk-c26g7-rabbit-mab/3333 Mucts/primary-antibodies/phospho-akt-thr308-antibody/9275 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9271 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9271 Aucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Aucts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/phospho-s6-ribosomal-protein-5g10-rabbit-mab/4695 aducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 aducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/12683 aducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3579 aducts/primary-antibodies/bash-acti-antibody/967 aducts/primary-antibodies/bash-acti-antibody/4967 aducts/primary-antibodies/anti-anbibi-igg/2729 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 aducts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/707
 https://www.cellsignal.com/prod https://www.cellsignal.com/produ https://www.cellsignal.com/produ https://www.cellsignal.com/produ https://www.eellsignal.com/produ https://www.eellsignal.com/produ https://www.eellsignal.com/produ https://www.eellsignal.com/produ https://www.abcam.com/alk-pd https://www.abcam.com/alk-pd https://www.abcam.com/bristor 	Aucts/primary-antibodies/n-myc-antibody/9405 Aucts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Bucts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 Bucts/primary-antibodies/alk-c26g7-rabbit-mab/4691 Bucts/primary-antibodies/phospho-akt-thr308-antibody/9275 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-set-ibosomal-protein-ser235-236-91b2-rabbit- abucts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 Bucts/primary-antibodies/phospho-st3-tyr705-antibody/9131 bucts/primary-antibodies/phospho-st3-tyr705-antibody/9131 bucts/primary-antibodies/bat3-79d7-rabbit-mab/4904 bucts/primary-antibodies/bactin-abcb1-d3h1q-rabbit-mab/12683 bucts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3579 bucts/primary-antibodies/ormal-rabbit-igg/2729 bucts/primary-antibodies/oromal-rabbit-igg/2729 bucts/primary-antibodies/normal-rabbit-igg/2729 bucts/primary-antibodies/normal-rabbit-igg/2729 bucts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 bucts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 bucts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 bucts/product/Anti-CTCF-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-CTCF-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-CTCF-Antibody,MM_NF-07-729 (US/en/product/Anti-CTCF-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-
 https://www.cellsignal.com/prod https://www.cellsignal.com/produ https://www.eellsignal.com/produ https://www.eellsignal.com/produ https://www.eellsignal.com/produ https://www.abcam.com/alk-pl https://www.abcam.com/histor https://www.novusbio.com/produ 	Aucts/primary-antibodies/n-myc-antibody/9405 Aucts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Bucts/primary-antibodies/alk-c26g7-rabbit-mab/3333 Aucts/primary-antibodies/alk-c26g7-rabbit-mab/4691 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9275 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- Bucts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 Bucts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 Bucts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 Bucts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 Bucts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 Bucts/primary-antibodies/bact-d3h1q-rabbit-mab/12683 Bucts/primary-antibodies/bactf-d1a7-xp-rabbit-mab/3579 Bucts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3579 Bucts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3417 Bucts/primary-antibodies/normal-rabbit-igg/2729 Bucts/secondary-antibodies/normal-rabbit-gg/2729 Bucts/secondary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 Bucts/primary-antibodies/anti-mouse-igg-hrp-linked-antibody/7076 Bucts/product/Anti-CTCF-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody-ADT3996.html a-ntibody-ep12204-ab187163.html he-h3-acetyl-k27-antibody-20b11_nbp2-52405
 https://www.cellsignal.com/prod https://www.eellsignal.com/prod https://www.eellsignal.com/prod https://www.eellsignal.com/prod https://www.active.com/ak-pd https://www.active.com/ak-pd https://www.active.com/prod https://www.active.com/prod https://www.active.com	Aucts/primary-antibodies/n-myc-dub2y-rabbit-mab/51705 Aucts/primary-antibodies/n-myc-dub2y-rabbit-mab/51705 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Aucts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 Aucts/primary-antibodies/alk-c26g7-rabbit-mab/4691 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9275 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9271 Aucts/primary-antibodies/phospho-akt-thr308-antibody/9271 Aucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Aucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Aucts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/phospho-s6-ribosomal-protein-seg10-rabbit-mab/2217 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/phospho-stat3-tyr705-antibody/9131 oducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/3579 oducts/primary-antibodies/cfcf-d1a7-xp-rabbit-mab/3579 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/primary-antibodies/normal-rabbit-igg/2729 oducts/product/Anti-Sox9-Antibody (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 (US/en/product/Anti-Sox9-Antibody,MM_NF-07-749 nospho-y1507-antibody-ab73996.html antibody-e012_04-ab187163.html he-h3-acetyl-k27-antibody-c0b11_nbp2-52405 tatalog/details/39851/boris-ctcfl-antibody-pab
 https://www.cellsignal.com/prod https://www.cellsignal.com/produ https://www.cellsignal.com/produ https://www.sethyl.com/produ https://www.sethyl.com/produ http://www.endmillipore.com/ http://www.abcam.com/ak-pl https://www.abcam.com/boris- https://www.abcam.com/boris- https://www.abcam.com/boris- https://www.abcam.com/boris- 	Aucts/primary-antibodies/n-myc-antibody/9405 Aucts/primary-antibodies/n-myc-d4b2y-rabbit-mab/51705 Aucts/primary-antibodies/cleaved-parp-asp214-antibody-human-specific/9541 Bucts/primary-antibodies/cleaved-caspase-3-asp175-antibody/9661 Bucts/primary-antibodies/alk-c26g7-rabbit-mab/4691 Bucts/primary-antibodies/phospho-akt-thr308-antibody/9275 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-akt-ser473-antibody/9271 Bucts/primary-antibodies/phospho-p44-42-mapk-erk1-2-thr202-tyr204-197g2-rabbit- aducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 Bucts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-s6-ribosomal-protein-ser235-236-91b2-rabbit-mab/4857 aducts/primary-antibodies/phospho-st3-tyr705-antibody/9131 aducts/primary-antibodies/phospho-st3-tyr705-antibody/9131 aducts/primary-antibodies/phospho-st3-tyr705-antibody/9131 aducts/primary-antibodies/phospho-st3-tyr705-antibody/9131 aducts/primary-antibodies/sox2-d6d9-xp-rabbit-mab/4579 aducts/primary-antibodies/cf-d1a7-xp-rabbit-mab/3579 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/primary-antibodies/normal-rabbit-igg/2729 aducts/primary-antibodies/normal-rabbit-igg/2729 /ust/A301-985A100/BRD4+Antibody /uS/en/product/Anti-CTCF-Antibody,MM_NF-07-729 /US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 /US/en/product/Anti-Sox9-Antibody,MM_NF-07-729 /US/en/product/Anti-Sox9-Antibody,MM_NF-07-749 nospho-y1507-antibody-ab73996.html antibody-ep12204-ab187163.html ne-h3-acetyl-k27-antibody-20b11_nbp2-52405

Validation

Eukaryotic cell lines

Policy information about <u>cell lines</u> Cell line source(s)	Human NB cell lines Kelly and SK-N-BE(2) and human Ewing sarcoma cell lines TC-32, TC-71 and CHLA-10 were obtained from the Children's Oncology Group cell line bank (Lubbock, TX, USA). Human NB cell line SK-N-SH and human embryonic kidney cell line HEK293T were obtained from the American Type Culture Collection (Manassas, VA, USA).
Authentication	All cell lines have been authenticated by short tandem repeat analysis.
Mycoplasma contamination	All cell lines tested negative for mycoplasma.
Commonly misidentified lines (See <u>ICLAC</u> register)	No commonly misidentified cell lines were used.

Animals and other organisms

Policy information about studies involving animals; ARRIVE guidelines recommended for reporting animal research

Laboratory animals	NU/NU (Crl:NU-Foxn1nu) mice (6-8 weeks old female) were purchased from Charles River Laboratories. NU/NU (CrTac:NCr-Foxn1nu) mice (6-8 weeks old female) were purchased from Taconic.
Wild animals	The study did not involve wild animals.
Field-collected samples	The study did not involve samples collected from the field.
Ethics oversight	All animal experiments were performed with approval from the Institutional Animal Care and Use Committee of the DFCI.

Note that full information on the approval of the study protocol must also be provided in the manuscript.

Human research participants

Policy information about studies involving human research participants

Population characteristics	Data involving human research participants pertain to immunohistochemistry of formalin-fixed paraffin-embedded tumor tissue slides, and were de-identified prior to analysis.
Recruitment	Data involving human research participants pertain to immunohistochemistry of formalin-fixed paraffin-embedded tumor tissue slides, and were de-identified prior to analysis.
Ethics oversight	All human tumor specimens (formalin-fixed paraffin-embedded slides) were obtained under an Institutional Review Board- approved protocol of the Dana-Farber/Boston Children's Cancer and Blood Disorders Center.

Note that full information on the approval of the study protocol must also be provided in the manuscript.

ChIP-seq

Data deposition

Confirm that both raw and final processed data have been deposited in a public database such as GEO.

Confirm that you have deposited or provided access to graph files (e.g. BED files) for the called peaks.

Data access links May remain private before publication.	https://www.ncbi.nlm.nih.gov/geo/query/acc.cgi?acc=GSE103084
Files in database submission	ChIPseq_Original_Submission_sens_res Raw files: Par_BORIS_1.fastq.gz Par_BORIS_2.fastq.gz Par_BRD4_1.fastq.gz Par_BRD4_2.fastq.gz Par_CTCF_1.fastq.gz Par_CTCF_2.fastq.gz Par_H3K27ac_1.fastq.gz Par_H3K27me3_1.fastq.gz Par_H3K27me3_2.fastq.gz Par_input_1.fastq.gz Par_input_2.fastq.gz Par_POI2_1.fastq.gz Par_POI2_2.fastq.gz

Res BORIS 1.fastq.gz Res_BORIS_2.fastq.gz Res_BRD4_1.fastq.gz Res_BRD4_2.fastq.gz Res_CTCF_1.fastq.gz Res_CTCF_2.fastq.gz Res_H3K27ac_1.fastq.gz Res_H3K27ac_2.fastq.gz Res_H3K27me3_1.fastq.gz Res_H3K27me3_2.fastq.gz Res_input_1.fastq.gz Res_input_2.fastq.gz Res_Pol2_1.fastq.gz Res_Pol2_2.fastq.gz Res_JQ1_BORIS_1.fastq.gz Res_JQ1_BORIS_2.fastq.gz Res_JQ1_BRD4_1.fastq.gz Res_JQ1_BRD4_2.fastq.gz Res_JQ1_CTCF_1.fastq.gz Hi2.fastq.gz Res_JQ1_H3K27ac_1.fastq.gz Res_JQ1_H3K27ac_2.fastq.gz Res_JQ1_H3K27me3_1.fastq.gz Res_JQ1_H3K27me3_2.fastq.gz Res_JQ1_input_1.fastq.gz Res_JQ1_input_2.fastq.gz Processed files. norm_Par_BORIS_merged.noNeg.sorted.bw norm_Par_BRD4_merged.noNeg.sorted.bw norm_Par_CTCF_merged.noNeg.sorted.bw norm_Par_H3K27ac_merged.noNeg.sorted.bw norm_Par_H3K27me3_merged.noNeg.sorted.bw norm_Par_Pol2_2.noNeg.sorted.bw norm_Res_BORIS_merged.noNeg.sorted.bw norm_Res_BRD4_merged.noNeg.sorted.bw norm Res CTCF merged.noNeg.sorted.bw norm_Res_H3K27ac_merged.noNeg.sorted.bw norm_Res_H3K27me3_merged.noNeg.sorted.bw norm_Res_Pol2_2.noNeg.sorted.bw norm_Res_JQ1_BORIS_merged.noNeg.sorted.bw norm_Res_JQ1_BRD4_merged.noNeg.sorted.bw norm_Res_JQ1_CTCF_merged.noNeg.sorted.bw norm_Res_JQ1_H3K27ac_merged.noNeg.sorted.bw norm_Res_JQ1_H3K27me3_merged.noNeg.sorted.bw Par_BORIS_vs_Par_input_peaks.narrowPeak Par_BRD4_vs_Par_input_peaks.narrowPeak Par_CTCF_vs_Par_input_peaks.narrowPeak Par_narrow_Pol2_vs_Par_input_peaks.narrowPeak Res_BORIS_vs_Res_input_peaks.narrowPeak Res_BRD4_vs_Par_input_peaks.narrowPeak Res_CTCF_vs_Res_input_peaks.narrowPeak Res_narrow_Pol2_vs_Res_input_peaks.narrowPeak Res_JQ1_BORIS_vs_Res_JQ1_input_peaks.narrowPeak Res_JQ1_BRD4_vs_Res_JQ1_input_peaks.narrowPeak Res_JQ1_CTCF_vs_Res_JQ1_input_peaks.narrowPeak Par_H3K27ac_vs_Par_input_peaks.broadPeak Par_H3K27me3_vs_Par_input_peaks.broadPeak Res_H3K27ac_vs_Res_input_peaks.broadPeak Res_H3K27me3_vs_Res_input_peaks.broadPeak Res JQ1 H3K27ac vs Res JQ1 input peaks.broadPeak Res_JQ1_H3K27me3_vs_Res_JQ1_input_peaks.broadPeak ChIPseq_MYCN Raw files: Kelly_input_MYCN_1.fastq.gz Kelly_MYCN_CS_1.fastq.gz Kelly_MYCN_CS_2.fastq.gz Res Kelly input MYCN CS 1.fastq.gz Res_Kelly_input_MYCN_CS_2.fastq.gz Res_Kelly_MYCN_CS_1.fastq.gz Res_Kelly_MYCN_CS_2.fastq.gz Processed files: Par_Kelly_MYCN.noNeg.sorted.bw Res_Kelly_MYCN.noNeg.sorted.bw Par_Kelly_MYCN_vs_Par_Kelly_MYCN_INPUT_peaks.narrowPeak Res_Kelly_MYCN_vs_Res_Kelly_MYCN_INPUT_peaks.narrowPeak ChIPseq KD Boris Raw files: shB4 Res Kelly BORIS 1.fastq.gz shB4_Res_Kelly_BORIS_2.fastq.gz shB4_Res_Kelly_CTCF_1.fastq.gz shB4_Res_Kelly_CTCF_2.fastq.gz shB4_Res_Kelly_H3K27ac_1.fastq.gz shB4_Res_Kelly_H3K27ac_2.fastq.gz shB4_Res_Kelly_INPUT_1.fastq.gz shB4_Res_Kelly_INPUT_2.fastq.gz shB4_Res_Kelly_INPUT_BORIS_1.fastq.gz shB4_Res_Kelly_INPUT_BORIS_2.fastq.gz shB4_Res_Kelly_SMC1_1.fastq.gz shB4_Res_Kelly_SMC1_2.fastq.gz shLUC_Res_Kelly_BORIS_1.fastq.gz shLUC_Res_Kelly_BORIS_2.fastq.gz shLUC_Res_Kelly_CTCF_1.fastq.gz shLUC_Res_Kelly_CTCF_2.fastq.gz shLUC_Res_Kelly_H3K27ac_1.fastq.gz shLUC_Res_Kelly_H3K27ac_2.fastq.gz shLUC_Res_Kelly_INPUT_1.fastq.gz shLUC_Res_Kelly_INPUT_2.fastq.gz shLUC_Res_Kelly_SMC1_1.fastq.gz shLUC_Res_Kelly_SMC1_2.fastq.gz Processed files: shB4_Res_Kelly_BORIS.noNeg.sorted.bw shB4_Res_Kelly_CTCF.noNeg.sorted.bw shB4_Res_Kelly_H3K27ac.noNeg.sorted.bw shB4_Res_Kelly_SMC1.noNeg.sorted.bw shLUC_Res_Kelly_BORIS.noNeg.sorted.bw shLUC_Res_Kelly_CTCF.noNeg.sorted.bw shLUC_Res_Kelly_H3K27ac.noNeg.sorted.bw shLUC_Res_Kelly_SMC1.noNeg.sorted.bw shB4_Res_Kelly_BORIS_vs_shB4_Res_Kelly_INPUT_peaks.narrowPeak shB4_Res_Kelly_CTCF_vs_shB4_Res_Kelly_INPUT_peaks.narrowPeak shB4 Res Kelly H3K27ac vs shB4 Res Kelly INPUT peaks.broadPeak shB4_Res_Kelly_SMC1_vs_shB4_Res_Kelly_INPUT_peaks.narrowPeak shLUC_Res_Kelly_BORIS_vs_shLUC_Res_Kelly_INPUT_peaks.narrowPeak shLUC_Res_Kelly_CTCF_vs_shLUC_Res_Kelly_INPUT_peaks.narrowPeak shLUC_Res_Kelly_H3K27ac_vs_shLUC_Res_Kelly_INPUT_peaks.broadPeak shLUC_Res_Kelly_SMC1_vs_shLUC_Res_Kelly_INPUT_peaks.narrowPeak ChIPseq_Ewing_Sarcoma Raw files: TC32_BORIS_1.fastq.gz TC32_BORIS_2.fastq.gz TC32_H3K27ac_1.fastq.gz TC32_H3K27ac_2.fastq.gz TC32_input_1.fastq.gz TC32_input_2.fastq.gz TC32_SMC1_1.fastq.gz TC32_SMC1_2.fastq.gz TC71 BORIS 1.fastq.gz TC71_BORIS_2.fastq.gz TC71_H3K27ac_1.fastq.gz TC71_H3K27ac_2.fastq.gz TC71_input_1.fastq.gz TC71_input_2.fastq.gz TC71_SMC1_1.fastq.gz TC71_SMC1_2.fastq.gz Processed files: Ewings_TC32_BORIS.noNeg.sorted.bw Ewings_TC32_H3K27ac.noNeg.sorted.bw Ewings_TC32_SMC1.noNeg.sorted.bw Ewings_TC71_BORIS.noNeg.sorted.bw Ewings_TC71_H3K27ac.noNeg.sorted.bw Ewings_TC71_SMC1.noNeg.sorted.bw Ewings TC32 BORIS vs David Ewings TC32 INPUT peaks.narrowPeak Ewings_TC32_H3K27ac_vs_David_Ewings_TC32_INPUT_peaks.broadPeak Ewings_TC32_SMC1_vs_David_Ewings_TC32_INPUT_peaks.narrowPeak Ewings_TC71_BORIS_vs_David_Ewings_TC71_INPUT_peaks.narrowPeak Ewings_TC71_H3K27ac_vs_David_Ewings_TC71_INPUT_peaks.broadPeak Ewings_TC71_SMC1_vs_David_Ewings_TC71_INPUT_peaks.narrowPeak Microarray_Original_Submission_sens_res_res_JQ1 Raw files: DD_1_PrimeView_par_Kelly_DMSO.CEL

DD_2_PrimeView_par_Kelly_DMSO.CEL DD_5_PrimeView_Res_Kelly_TAE.CEL DD_6_PrimeView_Res_Kelly_TAE.CEL DD_7_PrimeView_Res_Kelly_JQ1.CEL DD_8_PrimeView_Res_Kelly_JQ1.CEL Microarray_KD_Boris Raw files: Res_Kelly_shBORIS_1.CEL Res_Kelly_shBORIS_2.CEL Res_Kelly_shLUC_1.CEL Res_Kelly_shLUC_2.CEL 10X_sens_IR_Res Raw files: Sens.possorted_genome_bam.bam Intermediate_Res.possorted_genome_bam.bam Full_Res.possorted_genome_bam.bam Processed files: Sens_matrix.mtx Intermediate_Res_matrix.mtx Full_Res_matrix.mtx barcodes.tsv genes.tsv HiChIP_sens_res Raw files: KellyRes1_1.fastq.tar.gz KellyRes1_2.fastq.tar.gz KellyRes2_1.fastq.tar.gz KellyRes2_2.fastq.tar.gz KellyRes3_1.fastq.tar.gz KellyRes3_2.fastq.tar.gz KellyWTrep1_1.fastq.bz2 KellyWTrep1_2.fastq.bz2 KellyWTrep2_1.fastq.bz2 KellyWTrep2_2.fastq.bz2 KellyWTrep3_1.fastq.bz2 KellyWTrep3 2.fastq.bz2 Processed files: Kelly-Res-SMC1.bedpe Kelly-Res-SMC1-petcount.bedpe Kelly-SMC1.bedpe Kelly-SMC1-petcount.bedpe HiChIP_KD_Boris Raw files: Kelly_shBORIS1_1.fastq.tar.gz Kelly_shBORIS1_2.fastq.tar.gz Kelly_shBORIS2_1.fastq.tar.gz Kelly_shBORIS2_2.fastq.tar.gz Kelly_shBORIS3_1.fastq.tar.gz Kelly_shBORIS3_2.fastq.tar.gz Kelly_shLUC1_1.fastq.tar.gz Kelly_shLUC1_2.fastq.tar.gz Kelly_shLUC2_1.fastq.tar.gz Kelly_shLUC2_2.fastq.tar.gz Kelly_shLUC3_1.fastq.tar.gz Kelly_shLUC3_2.fastq.tar.gz Processed files: shBORIS-results.csv shLUC-results.csv

No longer applicable.

Genome browser session (e.g. <u>UCSC</u>)

Methodology


Replicates

Sequencing depth

For each mark assessed, two biological replicates were performed. All ChIP-seq data are derived from the analysis of both replicates.

ChIP-seq libraries were generated using the NEBNext Ultra II DNA Library Prep kit (E7645), following the manufacturer's instructions. Starting DNA material ranged from 2.5 to 10ng, and PCR amplification (8 to 10 cycles) was performed using NEBNext Multiplex Oligos for Illumina (E7335 and E7500) with distinct indices to allow for multiplexing of up to 12 samples to be run on the Illumina NextSeq 500 for 75 bases in single-read mode. Sample Total_reads Uniquely_mapped_reads Par_BORIS_1.fastq.gZ 52592735 39434783 Par_BORIS_2.fastq.gz 54031826 40040869

	Par_BRD4_1.fastq.gz 32913499 24145172
	Par_BRD4_2.fastq.gz 34899996 25050068
	Par_CTCF_1.fastq.gz 32357383 21742074
	Par_CTCF_2.fastq.gz 38048050 25570483
	Par_H3K27ac_1.fastq.gz 34652785 31481286
	Par_H3K27ac_2.fastq.gz 36652655 33723955
	Par_H3K27me3_1.fastq.gz 35375006 28384362
	Par_H3K27me3_2.fastq.gz 33791944 27655558
	Par_Pol2_1.fastq.gz 24373861 6139176
	Par_Pol2_2.fastq.gz 41084548 2685056 Kelly_MYCN_CS_1.fastq.gz 22544426 19447997
	Kelly_MYCN_CS_2.fastq.gz 43975115 38012488
	Par_input_1.fastq.gz 74638885 66630061
	Par_input_2.fastq.gz 40772400 36361135
	Kelly_input_MYCN_1.fastq.gz 43184072 34957425
	Res_BORIS_1.fastq.gz 43137736 32555678
	Res_BORIS_2.fastq.gz 47143635 35144692
	Res_BRD4_1.fastq.gz 32159271 24332339
	Res_BRD4_2.fastq.gz 36755542 30051197
	Res_CTCF_1.fastq.gz 40911369 28372298
	Res_CTCF_2.fastq.gz 36878692 25628082
	Res_H3K27ac_1.fastq.gz 33542051 30727706
	Res_H3K27ac_2.fastq.gz 37596270 34409319
	Res_H3K27me3_1.fastq.gz 35238601 29098617
	Res_H3K27me3_2.fastq.gz 35026938 28826263
	Res_Pol2_1.fastq.gz 35817388 6953374
	Res_Pol2_2.fastq.gz 48429363 2674818
	Res_Kelly_MYCN_CS_1.fastq.gz 18821385 14198653
	Res_Kelly_MYCN_CS_2.fastq.gz 30543884 24497611
	Res_Kelly_input_MYCN_CS_1.fastq.gz 77322913 68513199
	Res_Kelly_input_MYCN_CS_2.fastq.gz 53883106 47951287
	Res_input_1.fastq.gz 45072526 40107804 Res_input_2.fastq.gz 49535965 44608902
	Res_JQ1_BORIS_1.fastq.gz 81479030 61882578
	Res_JQ1_BORIS_2.fastq.gz 64057634 49039377
	Res_JQ1_BRD4_1.fastq.gz 32714245 25430455
	Res_JQ1_BRD4_2.fastq.gz 33031805 25335103
	Res_JQ1_CTCF_1.fastq.gz 33287924 24508642
	Res_JQ1_CTCF_2.fastq.gz 39674176 29154183
	Res_JQ1 H3K27ac 1.fastq.gz 40677082 37149902
	Res_JQ1_H3K27ac_2.fastq.gz 32078957 29293906
	Res_JQ1_H3K27me3_1.fastq.gz 36342163 30007069
	Res_JQ1_H3K27me3_2.fastq.gz 37930668 31338450
	Res_JQ1_input_1.fastq.gz 49313027 43993392
	Res_JQ1_input_2.fastq.gz 43284562 38433644
Antibodies	The following antibodies were used: MYCN (#51705, Cell Signaling Technology), BRD4 (A301-985A100, lot 6, Bethyl
	Laboratories), CTCF (#07-729, lot 2887267, Millipore), H3K27me3 (#07-449, lot 2972864, Millipore), H3K27ac (ab2729, lot
	GR3198866-1, Abcam), and BORIS (#39851, lot 18916002, Active Motif).
	All antibodies were previously validated by their manufacturers.
Peak calling parameters	Samples were aligned to the human genome (build hg19, GRCh37.75) with STAR (v2.5.1b modified) and the parameters "
reak cannig parameters	alignIntronMax 1alignEndsType EndToEndoutFilterMultimapNmax 1outFilterMismatchMax 5". Next, non-duplicate
	reads that mapped to the reference chromosomes were retained using Samtools (v1.3.1) and MarkDuplicates (v2.1.1) from
	Picard tools. Peaks were identified with MACS2 (2.1.1) for narrow peaks with the parameters "q 0.01call-summits" and
	for broad peaks with the parameters "broad-cutoff 0.01".
Data quality	1) Peaks overlapping regions with known artefact regions (http://mitra.stanford.edu/kundaje/akundaje/release/blacklists/)
	were blacklisted out.
	2) Antibody enrichment was assessed using the plotFingerprint command from deepTools (v2.2.4).
	3) Correlation of replicates was assessed with the deepTools command "multiBigwigSummary BED-file" using all bigwigs and
	identified peaks as input.
	4) broad peaks with 5-fold enrichment, 5% FDR and after blacklisting
	Par_H3K27ac_vs_Par_input_peaks.broadPeak.blacklisted 15626
	Par_H3K27me3_vs_Par_input_peaks.broadPeak.blacklisted 3946
	Res_H3K27ac_vs_Res_input_peaks.broadPeak.blacklisted 16725
	Res_H3K27me3_vs_Res_input_peaks.broadPeak.blacklisted 2293
	Res_JQ1_H3K27ac_vs_Res_JQ1_input_peaks.broadPeak.blacklisted 15942
	Res_JQ1_H3K27me3_vs_Res_JQ1_input_peaks.broadPeak.blacklisted 1819 5) narrow peaks with 5-fold enrichment, 5% FDR and after blacklisting
	Par_BORIS_vs_Par_input_peaks_blacklisted.narrowPeak 1055
	Par_BRD4_vs_Par_input_peaks_blacklisted.narrowPeak 52158
	Par_CTCF_vs_Par_input_peaks_blacklisted.narrowPeak 69566
	Par_narrow_Pol2_vs_Par_input_peaks_blacklisted.narrowPeak 19193
	Par_MYCN_vs_Par_input_peaks_blacklisted.narrowPeak 36047



Flow Cytometry

Plots

Confirm that:

The axis labels state the marker and fluorochrome used (e.g. CD4-FITC).

The axis scales are clearly visible. Include numbers along axes only for bottom left plot of group (a 'group' is an analysis of identical markers).

All plots are contour plots with outliers or pseudocolor plots.

A numerical value for number of cells or percentage (with statistics) is provided.

Methodology

Sample preparation	Flow cytometry was used for cell cycle analysis in both parental and resistant cells using propidium iodide (PI) staining. After fixation of 1x10e6 cells overnight at 4C with 80% ethanol, cells were resuspended in PBS supplemented with 0.1% Triton X-100, 25 mg/ml PI, and 0.2 mg/ml RNase A, and incubated for 45 min at 37C in the dark before analysis.
Instrument	Samples were run on a FACSCalibur flow cytometer (BD Biosciences).
Software	Data collection was done with BD CellQuest Pro software (BD Biosciences), and analysis was performed using FlowJo software (v10.0.5).
Cell population abundance	No FACS sorting was performed for this work.
Gating strategy	Cell debris as well as non-singlets were excluded by a gate based on FSC-H and SSC-H.

Tick this box to confirm that a figure exemplifying the gating strategy is provided in the Supplementary Information.